

# US580

LAPTOP ULTRASOUND  
DIAGNOSTIC SYSTEM



- \* Thin and light notebook design
- \* 33\*28.9\*7cm, 4.7KG
- \* 12.1 inch color LED display
- \* Rechargeable lithium battery
- \* Local zoom can be used at real time
- \* High cost performance always makes it a top seller

## ■ Various probes:

- 1).Convex probe
- 2).Linear probe
- 3).Transvaginal probe
- 4).Micro-convex probe
- 5).Rectal probe

Measurement:

General measurement package

GYN & OB measurement package

OB-1:

GS,BPD,CRL,FL,HL,TAD,LV,OFD,AC,HC,AIFI

OB-2:

FTA,TTD,APTD,THD,TCD,CI,EFW,Input LMP,Growth Curve

Andrology measurement package

Cardiac measurement package