

# USRKU 10 Veterinary Ultrasound Scanner

*(RING KING ULTRASOUND)*



USRKU 10 can be widely used in normal ultrasound diagnosis like abdomen and pregnancy.

Main market: big animals pregnancy check and follicular aspiration (Bovine, equine, ovine, canine, feline, goat, swine and llama).

## **Feature:**

- 1.1kg weigh. Smart and light.
- High-capacity lithium Battery ,more than 6 hours work.
- Sunshine cover, suitable for use in outdoor strong sunshine.
- Available for left/right hand alternative use, eliminating operate fatigue.
- Multiple use method: hand-held, chest hanging, back hanging , waist hanging, ect.
- Optional professional ultrasound goggle, convenient for outdoor use, meet different need.
- Have soft rubber coat, splash-proof, dust-proof.
- 2 kinds filter function, 7 kinds compress curve, 100 frames image permanent storage, 256 frames cine loop, Automatically calculate the Back fat and Lean percentage of swine.



## Tech Specification

- adopt 5.7" high resolution LED display, with High brightness, high contrast, wide viewing angle, suitable for outdoor use
- English language
- single probe socket, auto identify probe
- Body mark: 1.General 2. Bovine 3.Canine 4.Equine 5.Feline 6.Cardio 7. Swine, total 27 kinds body mark with probe director
- Display mode: B、2B、4B 、B/M、M
- Gain: 0~127dB
- Contrast: 27~90dB
- Gray: 256
- Focus: 1、2、3、4
- Magnification: 0.8、1.0、1.2、1.5、1.6、2.0、2.4、3.0
- Multi frequency: 4 different frequencies of each probe
- Depth: 8 steps, 240mm,
- Frame correlation: 0、1、2、3
- IP: 0、1、2、3



**Optional Configuration:**

- Probe:
  - 4.0MHz multi-frequency convex rectal Probe
  - 6.5MHz multi-frequency micro-convex probe
  - 3.5MHz multi-frequency convex probe
  - 6.5MHz multi-frequency trans-vaginal probe
  - 7.5MHz multi-frequency linear probe
  - 3.5MHz multi-frequency micro-convex cardiac probe
  - 3.5MHz loin muscle probe
- Thermal printer (P93W-S)
- Color laser printer (Specified models : HP Color LaserJet CP1515n、HP Color LaserJet CP1525n)
- Lithium Battery (5200mAH)
- Leather Case saddlebag