

WANROOEMED LONG TERM HOME CARE BAE509



Affordable comfort for residents

High comfort care bed with stepless motorized motion of the platform lifting system, the backrest and the foot rest. Suitable for use in nursing homes and in a homely environment.

It is used for families, hospitals etc, a popular product for EU market.

Features

1. The bed is designed for use with a standard mattress of 200 x 89 x 10 cm (not included)
2. 4-section platform, no-slip design.
3. actuators driven the motions of the platform, adjusted by the handset with single fault safety and locking device
4. **optional:**
A lifting (monkey) pole with triangular hand grip supports patients to easily lift into a sitting position
5. Wooden headboard and footboard
Full length **Wooden imitation aluminum alloy side rails** to protect people from falling off in raised position, and easy to get out when lowered.
6. 4 x 100 mm diameter individually lockable castors.
7. **optional: Emergency battery powered lowering in case of main power failure.**
8. Surface treatment: powder coated in Grey or Ivory.
9. Complied with IEC 60601-2-52-12-2009.

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663



BAE509 LONG TERM HOME CARE BED

- 10. All electrical components are IPX 4 rated.
- 11. Easy to assemble and convenient to transport.

When the side rail in the lowest position, small gap between the two siderails, more safety.
 And the upper side rail is only 7 cm above the bed surface



Mattress is optional



Extra Low type: 20-60cm height



Hand Controller with conductive switch in case of wrong operation



Wooden imitation aluminum alloy side rails

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663



BAE509 LONG TERM HOME CARE BED

Vertical Package

Bed frame and headboard/footboard	<u>siderails</u>

ordinary package

headboard/footboard	<u>siderails</u>	<u>Bed frame</u>

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663



BAE509 LONG TERM HOME CARE BED

Technical Parameters:

Platform dimension	2000X900mm				
Exterior dimension	2200X1050mm				
Max patient weight	180KGS				
Safe working load	205KGS				
Backrest adjustable range	0-70°				
Footrest adjustable range	0-40°				
Height adjustable range	390-800mm				
Trendelenburg positioning function	20°				
Reverse trendelenburg positioning function	20°				
Motor quantity and brand	4 sets, Chinese brand standard/ Dewart Motor optional				
Water-resistant	IPX4				
Distortion-free, highly stable side rails	Yes				
Package information	Vertical package: 20GP: 30pcs, 40GP: 65pcs, 40HQ:130pcs So if vertical package, the 40HQ is a best choice, as can load more goods!				
	item	Carton size	N.W.	G.W.	One bed/3 cartons
	Bed frame	960*420*1230mm	95	105	
	siderails	2010*155*110	4	6	
	Money pole-optional	1640*470*52mm	5	7	
	Ordinary package: 20GP: 44pcs, 40GP: 90pcs, 40HQ:100pcs Net weight: 105kgs, gross weight:117kgs				
	item	Carton size	N.W.	G.W.	One bed/4 cartons
	Headboard/footboard	1115*990*260mm	46	50	
	Bed frame	1115*990*297mm	46	50	
	siderails	2010*155*120	8	10	
Money pole-optional	1640*470*52mm	5	7		

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

BAE5091 LONG TERM HOME CARE BED

WANROOEMED ECONOMIC LONG TERM HOME CARE BAE5091



Comparing with BAE509 standard, only headboard/footboard different

Price is cheaper, quality is also good.



Wooden imitation aluminum alloy side rails

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

Electric Home Care Bed BAE5092 with Bed extension



Affordable comfort for residents

Guarantee: 3 years on Motor, 1 year on Bed

Features

1. pipe steel frame treated with electro-static spraying ,anti-aging and anti-rust.
2. Wooden board+Steel pipe head and footboard.
3. foldable steel side rails.
4. Foldable packing, can save storage space .
5. Five functions by hand controller
6. Richmat motor/Moteck Motor/ Dewart Motor
7. 4 castors with break
8. bed extension can be used for a taller use more safer and comfortable experience without loss the bed original feature.

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

Detailed Drawings:



Technical Parameters:

Product size	2100*960*310~710mm
Back rest lifting angle	0-70°
Knee break angle	0-50°
Reverse Trendelenburg	12°
Hi-lo adjustable	31~71cm

Standard accessories:

- 1) Four 5" castors all with brakes
- 2) 1 set of steel powder coated side rails
- 3) 1 set of detachable Head & Foot boards
- 4) 1 set of motors

Optional:

- 1) I.V. Pole
- 2) Mattress
- 3) Overbed table
- 4) Back-up battery
- 5) Bedside locker
- 6) lifting pole
- 7) extended leg section 18cm



no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

BAE05207 ELECTRIC HOME CARE BED

BAE05207 Electric Home Care Bed



Scope of application:

Product is suitable for the elderly, the sick, the patient mobility after surgery can not be moved or used under escort nursing staff, hospital, care center or family can be used.

Features:

Main frame made of high strength steel, surface powder spray, removable MDF headboard and footboard, removable guardrail two-stage, back lift, leg lift, the hand control two functions.

Specification:

Platform dimension	1900 X 850mm
Exterior dimension	2000X900mm
Safe working load	180KGS
Backrest adjustable range	0-70°
Footrest adjustable range	0-38°
Height adjustable range	430mm
Motor quantity and brand	2 sets, Chinese brand standard
Siderails height	400mm
Headboard height	460mm
Package information	20GP: 65pcs, 40GP: 130pcs, 40HQ:138pcs Net weight: 60kgs, gross weight:65kgs

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

Choose your surface

OPTIONAL 1:

FOAM HIGH-DENSITY SPONGE MATTRESS

SIZE: 2000*890MM, THICKNESS:100MM STANDARD

WATERPROOF COVER WITH ZIPPER



OPTIONAL 2:

MEDICAL AIR MATTRESS

1) **APP-B01 Anti-decubitus mattress**

APP-B01 is an anti-decubitus support system for medium risk patient

no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663



Mattress

Basic treatment therapeutic mattress supports is a good for pressure ulcer therapy. which includes an adjustable pump that allows for individual caregiver comfort settings. Heavy-duty, bubble-style pad is heat sealed by uniform tooling make pad more durable.

Key features:

bubbles on A-B lines alternating inflating provides superior therapeutic treatment

made of durable heavy , low toxic premium vinyl

Size: 1900*900mm

Cell height: 65mm

Material: Medical Grade PVC

Mattress thickness: 0.30mm

Number of cells:130pcs

Weight Support:120KG

Carton: 360*160*210mm

2300pcs/20GP,4600pcs/40GP,5700pcs/40HQ

2) APP-T01 Anti-decubitus mattress



APP-T01 is an anti-decubitus support system for medium risk patient

Key features:

modular air cells easy change and clean,
Convenience and Economic

vapor permeable ,water resistant and fire retardant cover

convert sealed tube durable, anti-friction
direction snaps make cover easy fix and
change

with ends flaps protect mattress slide on
bed

Size: 1900*900mm

Cell height: 90mm

Cells and Base Material: Nylon+PVC

Mattress thickness: 0.35mm

Number of cells:20 pcs

Max Weight Support:150KG

Carton: 370*345*168mm

1400pcs/20GP,2800pcs/4GP,

3300pcs/40HQ



no order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663