

ITEM		TECHNICAL REQUIREMENT
Size	Full Length	850mm
	Full Width	530mm
	Full Height	920mm
Castor Parameter		φ100mm Casters (with Diagonal Fixings)
Product Weight		16kg
Material	Bed Frame	High Quality Steel
	Sleeping Basin	Imported ABS
Range Of Motion	Angle	0~12°
Safe Working Load		Sleeping Basin: 15kg Basket: 5kg



Applicable medical and nursing institutions etc.

Medical baby bed

CARE AND LOVE ALWAYS WITH YOU



Frame

- ✦ The frame is made of high-quality steel-argon-arc welding and environmental protection coating, the safety performance of strong and corrosion resistant.




Caster

- ✦ Using $\phi 100\text{mm}$ double-sided casters, More stability and safety.



Basket

- ✦ Basket under the bed frame, more easy of Nursing.




Sleeping Basin

- ✦ Using high quality imported fully transparent ABS material.



Optional

- ✦ Mosquito net



Angle Adjusting Handle

- ✦ Using high quality gas spring to Stepless lift, Angle range 0-12°.



Mattress

- ✦ Fabric is waterproof and breathable, fire resistance, Flame retardant, antibacterial, non-toxic and so on.



Care And Love Always With You