

High Quality Drip Stand

BIV01 Adjustable IV stand(4hooks)

		
<p>4 legs metal base</p>	<p>5 legs plastic base</p>	<p>5 legs aluminum alloy base</p>

Pole:stainless steel 202#

Base:Powder-coated steel

Height adjustable, 4 hooks

Tube diameter: lower tube: 25mm, upper tube: 16 diameter, thickness: 1mm

110-180cm height adjustable

Plastic base diameter: 50cm

5 cm diameter castor with brakes.



BIV02 INFUSION PUMP RACK

1. Height adjustable easily
2. Stainless steel tube, with 4 IV hooks
3. original armrest help the patient adjust height according to the patient at any time and angle design principle of ergonomic comfort grip and lightweight operation.
4. Fast locking knob using Japanese technology, quick-locking and release, use one hand operation;
5. The base of the bracket with the universal self-Zohar pulley, flexible movement and braking.
6. With a plastic tray, applicable to all types of infusion pump equipment
7. Different kinds of bases optional



High Quality Drip Stand

BIV03 Adjustable Drip Stand

Pole:stainless steel 304

Rubber castors, 5 pcs

ABS base cover, base diameter: 450mm

Height adjustable:1260-1970mm, 2 hooks standard

Tube diameter: lower tube: 25mm, upper tube: 16 diameter

With Hand Rail, easy move the drip stand



BIV04 Drip Stand

Made in 304 stainless steel

Height adjustable, 4 hooks

Tube diameter: lower tube: 25mm, upper tube: 16 diameter, thickness: 1mm

110-180cm height adjustable

5 castors

