


Gynecologic and Obstetric
妇产计划生育



QUALITY
OUR WAY TO SUCCESS

Gynecologic and Obstetric



组织剪
TC METZENBAUM-FINE S-shaped scissors
ZC725RB
180 双弯 综合 镶片




组织剪
TC METZENBAUM scissors straight
ZC728RB
180 直 综合 镶片



子宫夹持钳
(子宫切除夹钳)
Hysterectomy Forceps
straight
ZG600RB
220 直 镀金圈
ZG607RB
250 直 镀金圈
ZG614RB
280 直 镀金圈

子宫夹持钳
(子宫切除夹钳)
Hysterectomy Forceps
slightly curved
ZG601RB
220×65×15 微弯 镀金圈
ZG608RB
250×78×15 微弯 镀金圈
ZG615RB
280×85×15 微弯 镀金圈

子宫夹持钳
(子宫切除夹钳)
Hysterectomy Forceps
curved
ZG602RB
220×62×20 中弯 镀金圈
ZG609RB
250×72×20 中弯 镀金圈
ZG616RB
280×81×20 中弯 镀金圈




子宫夹持钳
(子宫切除夹钳)
Hysterectomy Forceps
angled
ZG603RB
220×55×25 角弯 镀金圈
ZG610RB
250×65×25 角弯 镀金圈
ZG617RB
280×75×25 角弯 镀金圈

子宫夹持钳
(子宫切除夹钳)
Hysterectomy Forceps
angled
ZG604RB
220×90°×50 角弯 镀金圈
ZG611RB
250×90°×55 角弯 镀金圈
ZG618RB
280×90°×55 角弯 镀金圈

子宫夹持钳
(子宫切除夹钳)
Hysterectomy Forceps
S-shaped
ZG605RB
220×62×19 双弯 镀金圈
ZG612RB
250×72×23 双弯 镀金圈
ZG619RB
280×80×27 双弯 镀金圈


妇产计划生育
Gynecologic and Obstetric



子宫颈活体取样钳
Uterine Biopsy Forceps
3.5mm
ZG065R
180 长圆形
5.5mm
ZG066R
180, 圆头

子宫颈活体取样钳
Uterine Biopsy Forceps
3.5mm
ZG067RB
180 管式 长圆头
ZG069RB
230 管式 长圆头
5.5mm
ZG068RB
180 管式 圆头
ZG070RB
230 管式 圆头

妇产计划生育
Gynecologic and Obstetric




妇产计划生育
Gynecologic and Obstetric

妇产计划生育
Gynecologic and Obstetric




妇科钳
Duplay
tenaculum forceps
LS679R
280



子宫颈钳
Uterine forceps
curved
ZG176RB
180 弯 3×4齿
ZG177RB
220 弯 3×4齿
ZG178RB
250 弯 3×4齿



子宫颈钳
(宫颈抓取钳)
Uterine forceps
straight
ZG185R
250 3×4齿 直



子宫颈钳
Uterine forceps
straight
ZG179R
250 2×3齿 直
ZG179RN
250 2×3齿 直




子宫颈钳
Uterine forceps
side curved
ZG180R
250 2×3齿 侧弯
ZG180RN
250 2×3齿 侧弯



子宫息肉钳
WINTER Placenta
and Ovum Forceps
ZG303R
250×3 弯
ZG303RN
250×3 弯
ZG304R
250×6 弯
ZG304RN
250×6 弯



子宫息肉钳
WINTER Placenta
and Ovum Forceps
ZG300RB
250×6 直



子宫敷料钳
Bozeman Uterine
dressing forceps
ZG075R
250 无齿 双弯
ZG082R
250 有齿 双弯
ZG082RN
250 有齿 双弯



子宫敷料钳
Bozeman Uterine
dressing forceps, straight
ZG077RB
260 有齿 直



妇科组织钳
tissue forceps
ZE018R
200×5
ZG252RB
220×6
ZG253RB
240×6




止血钳
Artery forceps
ZH251RB
200 弯 竖齿



分离结扎钳
Ligature and
dissecting forceps
ZH290RB
210 角弯



妇产计划生育
Gynecologic and Obstetric



分离结扎钳
Ligature and
dissecting forceps
ZH292RB
230 弯90°

举宫钳
elevating forceps
ZG315RB
340 尖头

举宫钳
elevating forceps
ZG316RB
300 平头



妇产计划生育
Gynecologic and Obstetric



卵巢钳
Polypus forceps
straight
ZG317R
300 叠颚式 直

卵巢钳
Polypus forceps
curved
ZG318R
300 叠颚式 弯

卵巢钳
Polypus forceps
straight
ZG319RB
320 直
ZG320RB
320 弯



妇科分离钳
Dissecting Forceps
ZE053R
220×3 头圆弯R23×22



妇科分离钳
Dissecting Forceps
ZE054R
220×3 头圆弯R23×15



角弯阴道夹持钳
Vaginal holding forceps
angled
ZG956R
240×5 头圆弯R20×45



主韧带钳
Hysterectomy Forceps
ZG533R
220×4 微弯



泪滴状固定钳
Teardrop
Fixation Forceps
ZG469R
125



节育环放置钳
IUD applying
forceps
ZG020R
235

妇产计划生育
Gynecologic and Obstetric



主韧带钳
Hysterectomy Forceps
ZG518R
220×4 直



主韧带钳
Hysterectomy Forceps
ZG521R
220×4 头圆弯R23×22
ZG521RN
220×4 头圆弯R23×22



主韧带钳
Hysterectomy Forceps
ZG532R
210×5 头角弯22 竖齿1×2



宫内节育器取出钳
IUD removing forceps,curved
ZG003R
220 弯



宫内节育器取出钳
IUD removing forceps,straight
ZG016RG
200 直



附件钩钳
Uterine forceps
ZG183R
235 头圆弯R12×28

妇产计划生育
Gynecologic and Obstetric



产钳
Obstetric Forceps
ZG263RB
260 剖腹 单叶 折叠杆



产钳
Obstetric Forceps
ZG264R
300 单叶



产钳
Obstetric Forceps
ZG270R
350 低位



产钳
Obstetric Forceps
ZG265R
260 剖腹



产钳
Obstetric Forceps
ZG266RB
230 剖腹



产钳
Obstetric Forceps
ZG269R
400 转位



妇产计划生育
Gynecologic and Obstetric




妇产计划生育
Gynecologic and Obstetric



子宫动脉夹持钳
Uterine Tenaculum Forceps
curved
ZG181R
220 弯



输精管皮外固定钳 (固定钳)
Vasectomy fixation forceps
ZJ196R
140×2.2
ZG470R
140×3
ZG472R
140×3.5



阑尾肠钳
Babcock tissue forceps
ZE032RB
a
160×8.5 直
ZE033RB
200×11 直
ZE034RB
240×11 直



剖腹产切口钳
Caesarian Forceps
ZG271R
200 直 开口2
ZG272R
200 弯 开口2



输尿管夹持钳
tissue forceps
ZG471R
200 头带三角孔
ZG471RN
200 头带三角孔



输精管分离钳
Dissecting Forceps
ZJ197R
125 无齿



组织镊
tissue forceps
ZD076RB
230 直 带齿



环形输卵管镊
NEUBUESER Tubal
Seizing FORCEPS
ZD801R
200×13
ZD801RN
200×13



环形输卵管镊
NEUBUESER Tubal
Seizing FORCEPS
ZD818RN
160×8



环形输卵管镊
NEUBUESER Tubal
Seizing FORCEPS
ZD820R
200×Φ6
ZD821R
200×Φ8
ZD822R
200×Φ10





子宫拉钩
Vaginal Retractor
ZT613R
a b
250×40×114 平



子宫拉钩
Vaginal Retractor
ZT623RB
180×13×56 凹 扁柄



子宫拉钩
Vaginal Retractor
ZT624RB
180×10×60 平 扁柄





阴道拉钩
Vaginal Retractor

- ZG537RB** 220×14×115 平板式
- ZG538RB** 220×20×115 平板式
- ZG594RB** 220×24×115 平板式
- ZG539RB** 220×26×115 平板式
- ZG540RB** 220×32×115 平板式

阴道拉钩
Vaginal Retractor

- ZG595RB** 300×25×100 凹柄式
- ZG596RB** 300×30×100 凹柄式
- ZG565RB** 310×25×100 凹柄式
- ZG568RB** 320×40×100 凹柄式
- ZG597RB** 330×25×130 凹柄式
- ZG566RB** 330×30×100 凹柄式
- ZG598RB** 330×30×130 凹柄式
- ZG599RB** 350×25×150 凹柄式
- ZG567RB** 350×35×100 凹柄式



阴道拉钩
Vaginal Retractor

- ZG187R** $\begin{matrix} a & b \\ 285 \times 30 \times 85 \end{matrix}$ 凹柄式
- ZG188R** 300×33×100 凹柄式
- ZG189R** 315×36×115 凹柄式



阴道拉钩
Vaginal Retractor
ZG577R
240×40×100
重锤 1kg



阴道拉钩
Vaginal Retractor
ZG570RB
240×40×105
凹柄式 1.5kg重锤



宫颈压板
Cervical clamp
ZG222R
250×32 弯



阴道拉钩
Vaginal Retractor

- ZG571R** $\begin{matrix} a & b \\ 180 \times 30 \times 100 \end{matrix}$
- ZG571RN** 180×30×100
- ZG572R** 180×32×115
- ZG572RN** 180×32×115
- ZG573R** 180×36×126
- ZG573RN** 180×36×126



阴道压板
Vaginal Retractor

- ZG584R** 260×30×25° 凹柄式
- ZG586R** 220×30×45° 凹柄式



阴道压板
Vaginal Retractor

- ZG223RB** 250×30×35° 平板式



妇产计划生育
Gynecologic and Obstetric



妇产计划生育
Gynecologic and Obstetric

深部拉钩
Retractor
ZT654RB
260×35×90 伞柄



深部拉钩
Retractor
ZT658RB
260×50×120 扁柄
ZT659RB
260×39×150 扁柄
ZT653RB
260×57×150 扁柄
ZT660RB
250×35×91 扁柄
ZT661RB
260×50×180 扁柄



深部拉钩
Retractor
ZT655RB
290×45×130 伞柄
ZT657RB
290×60×130 伞柄



深部拉钩
Retractor
ZT662RB
320×25×160 板式
ZT649RB
320×50×160 板式
ZT665RB
360×25×180 板式
ZT651RB
360×50×180 板式




深部拉钩
Deaver retractor
ZT003RB
S形 330×25×100
ZT004RB
S形 300×38×100
ZT005RB
S形 300×50×100



耻骨联合上拉钩
Retractor
ZG581R
240×80×50



膀胱拉钩
Retractor
ZT496RB
240×75×45



拉钩
(输尿管拉钩)
Retractor
ZT162RB
静脉 200×16
ZT163RB
静脉 200×18


拉钩 (宫颈定位钩)
fistula hook
ZG642RB
220 尖头

结扎缝合引线器
(带线钩)
suture guide
left curved, blunt
ZG641RB
240 左弯 钝头


结扎缝合引线器
(带线钩)
suture guide,
right curved, blunt
ZG640RB
240 右弯 钝头

子宫刮匙
Uterine Curette

a



锐




- a
- ZG350RB** 300×5 锐 八角柄 柔性
 - ZG351RB** 300×7 锐 八角柄 柔性
 - ZG352RB** 300×8 锐 八角柄 柔性
 - ZG354RB** 300×10 锐 八角柄 柔性
 - ZG356RB** 300×12 锐 八角柄 柔性
 - ZG357RB** 300×14 锐 八角柄 柔性
 - ZG391RB** 300×15 锐 八角柄 柔性
 - ZG358RB** 300×16 锐 八角柄 柔性
 - ZG369RB** 300×18 锐 八角柄 柔性
 - ZG359RB** 300×20 锐 八角柄 柔性
 - ZG370RB** 300×22 锐 八角柄 柔性
 - ZG371RB** 300×24 锐 八角柄 柔性

钝




- a
- ZG392RB** 300×7.5 钝 八角柄 柔性
 - ZG360RB** 300×9 钝 八角柄 柔性
 - ZG361RB** 300×10 钝 八角柄 柔性
 - ZG363RB** 300×12 钝 八角柄 柔性
 - ZG364RB** 300×14 钝 八角柄 柔性
 - ZG393RB** 300×15 钝 八角柄 柔性
 - ZG365RB** 300×16 钝 八角柄 柔性
 - ZG372RB** 300×18 钝 八角柄 柔性
 - ZG366RB** 300×20 钝 八角柄 柔性
 - ZG394RB** 300×21.5 钝 八角柄 柔性
 - ZG395RB** 300×23.5 钝 八角柄 柔性

锐



- a
- ZG379RB** 255×7.5 锐 柔性 扁柄
 - ZG380RB** 255×8.5 锐 柔性 扁柄
 - ZG381RB** 255×10.5 锐 柔性 扁柄
 - ZG382RB** 255×12 锐 柔性 扁柄
 - ZG383RB** 255×14 锐 柔性 扁柄
 - ZG384RB** 255×15 锐 柔性 扁柄

钝



- a
- ZG385RB** 255×7.5 钝 柔性 扁柄
 - ZG386RB** 255×8.5 钝 柔性 扁柄
 - ZG387RB** 255×10.5 钝 柔性 扁柄
 - ZG388RB** 255×12 钝 柔性 扁柄
 - ZG389RB** 255×14 钝 柔性 扁柄
 - ZG390RB** 255×15 钝 柔性 扁柄






子宫刮匙
Uterine Curette

锐




钝



- | | |
|-------------------------------------|-------------------------------------|
| ZG337R 280×4 (4#) 锐 六方柄 柔性 | ZG401R 280×5 (5#) 钝 六方柄 柔性 |
| ZG338R 280×5 (5#) 锐 六方柄 柔性 | ZG402R 280×6 (6#) 钝 六方柄 柔性 |
| ZG397R 280×6 (6#) 锐 六方柄 柔性 | ZG403R 280×7 (7#) 钝 六方柄 柔性 |
| ZG339R 280×7 (7#) 锐 六方柄 柔性 | ZG404R 280×9 (9#) 钝 六方柄 柔性 |
| ZG398R 280×8 (8#) 锐 六方柄 柔性 | ZG405R 280×10 (10#) 钝 六方柄 柔性 |
| ZG340R 280×9 (9#) 锐 六方柄 柔性 | ZG348R 280×11 (11#) 钝 六方柄 柔性 |
| ZG399R 280×10 (10#) 锐 六方柄 柔性 | ZG406R 280×12 (12#) 钝 六方柄 柔性 |
| ZG341R 280×11 (11#) 锐 六方柄 柔性 | ZG408R 280×14 (14#) 钝 六方柄 柔性 |
| ZG400R 280×12 (12#) 锐 六方柄 柔性 | ZG407R 280×15 (15#) 钝 六方柄 柔性 |
| ZG342R 280×14 (14#) 锐 六方柄 柔性 | ZG409R 280×16 (16#) 钝 六方柄 柔性 |
| ZG343R 280×16 (16#) 锐 六方柄 柔性 | ZG410R 280×18 (18#) 钝 六方柄 柔性 |
| ZG344R 280×18 (18#) 锐 六方柄 柔性 | ZG349R 280×20 (20#) 钝 六方柄 柔性 |
| ZG345R 280×20 (20#) 锐 六方柄 柔性 | ZG433R 280×22 (22#) 钝 六方柄 柔性 |
| ZG346R 280×22 (22#) 锐 六方柄 柔性 | ZG432R 280×24 (24#) 钝 六方柄 柔性 |
| ZG347R 280×24 (24#) 锐 六方柄 柔性 | |

- | |
|-------------------------------------|
| ZG401R 280×5 (5#) 钝 六方柄 柔性 |
| ZG402R 280×6 (6#) 钝 六方柄 柔性 |
| ZG403R 280×7 (7#) 钝 六方柄 柔性 |
| ZG404R 280×9 (9#) 钝 六方柄 柔性 |
| ZG405R 280×10 (10#) 钝 六方柄 柔性 |
| ZG348R 280×11 (11#) 钝 六方柄 柔性 |
| ZG406R 280×12 (12#) 钝 六方柄 柔性 |
| ZG408R 280×14 (14#) 钝 六方柄 柔性 |
| ZG407R 280×15 (15#) 钝 六方柄 柔性 |
| ZG409R 280×16 (16#) 钝 六方柄 柔性 |
| ZG410R 280×18 (18#) 钝 六方柄 柔性 |
| ZG349R 280×20 (20#) 钝 六方柄 柔性 |
| ZG433R 280×22 (22#) 钝 六方柄 柔性 |
| ZG432R 280×24 (24#) 钝 六方柄 柔性 |




子宫刮匙 (宫颈刮匙)
Uterine Curette


ZG480RB
300

子宫刮匙
Uterine Curette

锐



钝



- | | |
|---------------------------------|-------------------------------|
| ZG411R 280×6 (6#) 锐 柔性 | ZG417R 280×5 (5#) 钝 柔性 |
| ZG412R 280×8 (8#) 锐 柔性 | ZG418R 280×6 (6#) 钝 柔性 |
| ZG413R 280×10 (10#) 锐 柔性 | ZG419R 280×7 (7#) 钝 柔性 |

- | |
|---------------------------------|
| ZG417R 280×5 (5#) 钝 柔性 |
| ZG418R 280×6 (6#) 钝 柔性 |
| ZG419R 280×7 (7#) 钝 柔性 |
| ZG420R 280×9 (9#) 钝 柔性 |
| ZG421R 280×10 (10#) 钝 柔性 |
| ZG422R 280×12 (12#) 钝 柔性 |
| ZG423R 280×15 (15#) 钝 柔性 |
| ZG424R 280×14 (14#) 钝 柔性 |
| ZG425R 280×16 (16#) 钝 柔性 |
| ZG426R 280×18 (18#) 钝 柔性 |



子宫探针
Uterine probe
ZG011R
280 直

- 子宫探针
Uterine probe
- ZG013RB** 320×Φ3 弯
 - ZG014RB** 320×Φ4 弯
 - ZG015RB** 320×Φ5 弯



会阴牵开器
Retractor
ZG492RB
180 夹式



会阴牵开器
Gelpi retractor
ZG490RB
170 短钩式



会阴牵开器
Gelpi retractor
ZG491RB
180 长钩式



双翼阴道手术扩张器
GRAVES Vaginal Speculum
ZG042R
92×34

双翼阴道检查扩张器
CUSCO Vaginal Speculum
ZG043R
90×22
ZG044R
100×32



拉钩
Retractor
ZV530R
200×100×105, 大号



拉钩
Retractor
ZV531R
200×90×85, 小号



腹壁牵开器
Retractor
ZV517R
280, 开口300, 三翼, 手摇式



腹壁牵开器
Retractor
ZV529R
最大开口 210, 钳式, 三翼, 固定



腹壁牵开器
Retractor
ZV520R
II型 (260×190)
钳式, 三翼, 固定



腹壁牵开器
Retractor
ZV579R
双圆梗, 175×202, 双翼, 固定

妇产计划生育
Gynecologic and Obstetric

妇产计划生育
Gynecologic and Obstetric




流产吸引管
Abortion
suction tube
ZG705C, 5#
ZG706C, 6#
ZG707C, 7#
ZG708C, 8#
ZG705RB, 5#
ZG706RB, 6#
ZG707RB, 7#
ZG708RB, 8#
ZG709RB, 9#
ZG713RB, 10#

腹壁牵开器
Retractor
ZV540R
175×202, 三翼, 固定



子宫颈扩张器
Uterine dilator, round tip

- | | | |
|--------------------|---------------------|--------------------|
| ZG273RB
φ3 圆头 | ZG284RB
φ8.5 圆头 | ZG295RB
φ16 圆头 |
| ZG274RB
φ3.5 圆头 | ZG285RB
φ9 圆头 | ZG296RB
φ17 圆头 |
| ZG275RB
φ4 圆头 | ZG286RB
φ9.5 圆头 | ZG297RB
φ18, 圆头 |
| ZG276RB
φ4.5 圆头 | ZG287RB
φ10 圆头 | ZG298RB
φ19 圆头 |
| ZG277RB
φ5 圆头 | ZG288RB
φ10.5 圆头 | ZG305RB
φ20 圆头 |
| ZG278RB
φ5.5 圆头 | ZG289RB
φ11 圆头 | ZG306RB
φ21 圆头 |
| ZG279RB
φ6 圆头 | ZG290RB
φ11.5 圆头 | ZG307RB
φ22 圆头 |
| ZG280RB
φ6.5 圆头 | ZG291RB
φ12 圆头 | ZG308RB
φ23 圆头 |
| ZG281RB
φ7 圆头 | ZG292RB
φ13 圆头 | ZG309RB
φ24 圆头 |
| ZG282RB
φ7.5 圆头 | ZG293RB
φ14 圆头 | |
| ZG283RB
φ8 圆头 | ZG294RB
φ15 圆头 | |



子宫颈扩张器
Uterine dilator

- | | |
|--------------------|---------------------|
| ZG090R
3# 尖圆头 | ZG113R
9# 尖圆头 |
| ZG102R
3.5# 尖圆头 | ZG117R
9.5# 尖圆头 |
| ZG104R
4# 尖圆头 | ZG114R
10# 尖圆头 |
| ZG103R
4.5# 尖圆头 | ZG118R
10.5# 尖圆头 |
| ZG105R
5# 尖圆头 | ZG115R
11# 尖圆头 |
| ZG106R
5.5# 尖圆头 | ZG092R
11.5# 尖圆头 |
| ZG107R
6# 尖圆头 | ZG093R
12# 尖圆头 |
| ZG108R
6.5# 尖圆头 | ZG094R
13# 尖圆头 |
| ZG109R
7# 尖圆头 | ZG095R
14# 尖圆头 |
| ZG110R
7.5# 尖圆头 | ZG096R
15# 尖圆头 |
| ZG111R
8# 尖圆头 | ZG097R
16# 尖圆头 |
| ZG112R
8.5# 尖圆头 | |



骨盆测量计
DE LEE Pelvimeter
ZG005R
200 测量范围10-120



骨盆测量计
(耻骨角度测量器)
Pubis Angle Meter
ZG004R
140 耻骨角度



骨盆测量计
DE LEE Pelvimeter
ZG007RB
330, 内/外径



妇产计划生育
Gynecologic and Obstetric



宫内节育器放置叉
IUD applier
ZG001C
280

宫内节育器取出钩
IUD removing hook
ZG002C
280

输卵管剥离器
dissector
ZG473R
220

输精管拉钩
hook
ZJ198R
135

输卵管提取钩
retractor
ZG475R
280×7.5×10



输卵管通液器
tubal irrigation catheter
ZG477C
320

输卵管提取板
tubal retractor
ZG479RB
150

肌瘤剥离器
Myoma dissector
ZG583R
320×12×12

破膜钩
Beacham Amniotome
ZG588R
250



妇产计划生育
Gynecologic and Obstetric



妇产计划生育
Gynecologic and Obstetric



结扎导板
Retractor
ZG958R
150 50×12 35×12

指板拉钩
Retractor
ZG959R
210 60×13 37×13

压肠板
Abdominal Spatula
ZT758R
300×35×45
ZT759R
350×35×45

压肠板
Abdominal Spatula
ZT760R
260×28, 带柄