

support you care you



WANROOEMED

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY

AOT400



AOT400A



AOT400C

support you care you

AOT400

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY



Introduction

Endearing appearance, various colors of mattress available., optional pedal switch electrically and easily operate.

It has overcome the slowness matter of similar products.

Main Technical Data

Parameters may have reasonable error.

AOT400A	AOT400C
Length and width for tabletop : 1820/520mm	Length and width for tabletop : 1780/600mm
Elevation(down/up) : 680/1090mm	Elevation(down/up) : 680/880mm
Trendelenburg/anti-trendelenburg : 14°/17°	Trendelenburg/anti-trendelenburg : 8°/20°
Back plate(up/down):60°/10°	Back plate(up/down):45°/7°
Power : 220V , 50Hz,1.0kW	Power : 220V , 50Hz,1.0Kw

support you care you



WANROOEMED

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY

AOT400B、AOT400BM

AOT400B



AOT400BM



support you care you

AOT400B、AOT400BM

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL
SURGERY



Introduction

400B adopted both remote control and pedal switch, diversified operation.

The cover for base and column is made of 304 stainless steel, good performance on anti-rust and easy to clean. 400BM with endearing appearance, seamless mattress, mechanically hydraulic control. Backplate is adjusted by germany made gas spring.

Main Technical Data

Parameters may have reasonable error.

AOT400B	AOT400BM
Length and width for tabletop : 1780/600mm	Length and width for tabletop : 1780/600mm
Elevation(down/up) : 680/880mm	Elevation(down/up) : 680/880mm
Trendelenburg/anti-trendelenburg : 8°/20°	
Back plate(up/down):45°/7°	Back plate(up/down):45°/10°
Power : 220V,50Hz,1.0kW	