

# 卓越而生，蓄势待发， ——感受万容轮椅之美学

Pursuing excellent and gaining momentum,  
to feel the aesthetics of Wanroo Wheelchair.



精妙的创意，精准人性化设计，精密轮椅技术科技智能核心、精湛工艺实现尊贵与舒适。  
优越生活，坐享健康，尊享未来。

Delicacy and full of originality design, accurate and human-based wheelchair technique.  
Use high-tech and intelligent crystal, creating a series of high-performance wheelchair product, bring people noble and comfortable enjoyment. Make a superior life for you. Let you appreciate the health and get beautiful feelings.

## 目录 Contents

钢质手动轮椅系列 Steel Manual Wheelchair Series	01~11
铝合金轮椅系列 Aluminum Wheelchair Series	11~18
电动轮椅系列 Electric Wheelchair Series	19
三轮车系列 3-wheel Rollator Series	20
手/拐/腋杖系列 Walking Stick / Crutch/ Armpit Crutch Series	21~22
助行器系列 Walking Aid Series	23
座厕椅系列 Commode Chair Series	24
护理床系列 Care Bed Series	25~26



## 钢质手动轮椅系列 Steel Manual Wheelchair Series



**BWHM-1B** 普通轮椅 General wheelchair

镀铬车架，单工架结构，固定扶手和搁脚，牛津布座靠垫，8"实芯前轮，24"免充气后轮，高强度塑料脚踏板。

Chromed frame, single frame structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" free pneumatic rear wheel, high strength plastic pedals.

64	50	39	43	29	45	104	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	26.5	17	96x26x92			

**BWHM-1B1** 普通轮椅 General wheelchair



镀铬车架，双交叉杆结构，固定扶手和搁脚，硬座，人造革座靠垫，8"实芯前轮，24"充气后轮，高强度塑料脚踏板。

Chromed frame, double cross bar structure, fixed armrest and footrests, hard seat, imitation-leather seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals.

64	50	39	43	23	45	104	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	20	92x26x95			

**BWHM-1B2** 普通轮椅 General wheelchair



**BWHM-1A** 多功能型轮椅 Multifunctional wheelchair

镀铬车架，单交叉杆结构，固定扶手，可拆卸搁脚，牛津布座靠垫，8"实芯前轮，24"免充气后轮。

Chromed frame, single cross bar structure, fixed armrest, detachable footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" free pneumatic rear wheel.

64	50	39	43	29	45	104	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	26.5	17	78x24x90			



喷塑车架，双交叉杆结构，固定扶手和搁脚，牛津布座靠垫，8"实芯前轮，24"充气后轮，高强度塑料脚踏板。

Plastic-spray frame, double cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals.

64	49	40	43	24	47	104	90
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	17	96x26x92			



## 钢质手动轮椅系列 Steel Manual Wheelchair Series



**BWHM-1B201** 普通轮椅 General wheelchair

喷塑车架，双交叉杆结构，固定扶手和搁脚，牛津布座靠垫，8" 实芯前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Chromed frame, double cross bar structure, fixed armrest and footrests, hard seat, imitation-leather seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

64	49	40	43	24	47	104	90
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	17	96x26x92			

**BWHM-1B14** 普通轮椅 General wheelchair



喷塑车架，单交叉杆结构，固定扶手和搁脚，牛津布座靠垫，8" 实芯前轮，24" 充气后轮，高强度塑料脚踏板。

Plastic-spray frame, single cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals.

64	50	39	43	29	45	104	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	26.5	17	96x26x92			

**BWHM-1B14** 普通轮椅 General wheelchair



**BWHM-1B11** (镀铬车架)

**BWHM-1B12** 普通轮椅 General wheelchair

喷塑车架，双交叉杆结构，固定扶手和搁脚，硬座，牛津布座靠垫，8" 实芯前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Plastic-spray frame, double cross bar structure, fixed armrest and footrests, hard seat, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

64	50	39	43	23	45	104	90
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	20	92x26x95			



喷塑车架，单交叉杆结构，固定扶手和搁脚，牛津布座靠垫，8" 实芯前轮，24" 充气后轮，带后刹车，高强度塑料脚踏板，附安全带。

Plastic-spray frame, single cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, with the rear brake, high strength plastic pedals, with safety belt.

64	50	39	43	29	45	104	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	26.5	17	96x26x92			



## 钢质手动轮椅系列 Steel Manual Wheelchair Series

### BWHM-1B29 加宽型轮椅 Widened wheelchair



喷塑车架，双交叉杆结构，固定扶手和搁脚，硬座，牛津布座靠垫，8" 实芯前轮，24" 充气后轮，高强度塑料脚踏板，座垫宽度为51cm，可载重150KG。

Plastic-spray frame, double cross bar structure, fixed armrest and footrests, hard seat, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals, seat width is 51cm, which can load 150KG.

70.5	50	39	43	23	51	104	90
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	20	97x28x99			

### BWHM-1A7 背部可调角度 Backside angle adjustable



镀铬车架，双交叉杆结构，高靠背，半躺，硬座，可拆卸搁脚，8" 加厚前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Chromed frame, double cross bar structure, high back, half lying, hard seat, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

66	52.5	41	46	20.5	45	114	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	24	92x26x95			

### BWHM-1A8 背部可调角度 Backside angle adjustable



### BWHM-1A5 背部可调角度 Backside angle adjustable

镀铬车架，双交叉杆结构，高靠背，半躺，可拆卸高抬搁脚，8" 加厚前轮，24" 充气后轮，附安全带。

Chromed frame, double cross bar structure, high back, half lying, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, with safety belt.

66	52.5	41	46	20.5	45	114	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	24	92x26x95			



镀铬车架，双交叉杆结构，高靠背，半躺，硬座，可拆卸高抬搁脚，8" 加厚前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Chromed frame, double cross bar structure, high back, half lying, hard seat, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

66	52.5	41	46	20.5	45	114	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	24	92x26x95			



## 钢质手动轮椅系列 Steel Manual Wheelchair Series



### BWHM-1A19 多功能型轮椅 Multifunctional wheelchair

烤漆车架，双交叉杆结构，翻转扶手，可拆卸搁脚，牛津布座靠垫，8" 铝合金前轮，24" 充气后轮，附安全带。

Stoving varnish frame, double cross bar structure, turnover armrest, detachable footrests, oxford cloth seat and backrest cushion, 8" aluminum alloy front wheel, 24" pneumatic rear wheel, with safety belt.

64	51	41	42	22.5	45	103	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	28	19	78x26x91			



### BWHM-1B502 (木侧板)

多功能型轮椅 Multifunctional wheelchair

喷漆车架，木质座厕板，附无孔PVC座垫，圆形灰色有盖深马桶，8" 实芯前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Plastic-spray frame, wood commode board, with imperforated PVC seat, rounded gray covered deep toilet, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

64	48	42	43	24	43	104	86
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	19	92x26x90			

### BWHM-1B5 (泡沫座侧垫) 座厕轮椅 Commode wheelchair



### BWHM-1B5 (木侧板)

座厕轮椅 Commode wheelchair

镀铬车架，木质座厕板，附无孔PVC座垫，圆形灰色有盖深马桶，8" 实芯前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Chromed frame, wood commode board, with imperforated PVC seat, rounded gray covered deep toilet, 8" solid front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

64	48	42	43	24	43	104	86
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	19	92x26x90			



镀铬车架，固定扶手和搁脚，人造革座靠垫，可抽拉座厕垫，8" 实芯前轮，24" 充气后轮。

Chromed frame, fixed armrest and footrests, imitation-leather seat and backrest cushion, pulling commode cushion, 8" solid front wheel, 24" pneumatic rear wheel.

64	48	42	43	24	43	104	86
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	19	94x28x91			



## 钢质手动轮椅系列 Steel Manual Wheelchair Series



### BWHM-1B502 (泡沫座厕垫)

座厕轮椅 Commode wheelchair

喷塑车架，固定扶手和搁脚，人造革座靠垫，可抽拉座厕垫，8" 实芯前轮，24" 充气后轮。

Plastic-spray frame, fixed armrest and footrests, imitation-leather seat and backrest cushion, pulling commode cushion, 8" solid front wheel, 24" pneumatic rear wheel.

64	48	42	43	24	43	104	86
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	19	94x28x91			



### BWHM-1A8-2 多功能型轮椅 Multifunctional wheelchair

喷塑车架，双交叉杆结构，高靠背，半躺，改进型可抽拉座厕垫，防后翻装置，后把手连刹车，靠背可按需要调节成不同的角度，可拆卸高抬搁脚，8" 加厚前轮，24" 充气后轮，附安全带。

Plastic-spray frame, double cross bar structure, high back, half lying, improved pulling commode cushion, anti overturn device, the rear handle link the brakes, backrest can adjust the different angle, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, with safety belt.

66	52.5	44.5	45	20.5	43	118	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	26	92x26x95			

### BWHM-1A902 座厕轮椅 Commode wheelchair



### BWHM-1A8-1 座厕轮椅 Commode wheelchair

镀铬车架，双交叉杆结构，高靠背，半躺，改进型可抽拉座厕垫，防后翻装置，后把手连刹车，靠背可按需要调节成不同的角度，可拆卸高抬搁脚，8" 加厚前轮，24" 充气后轮，附安全带。

Chromed frame, double cross bar structure, high back, half lying, improved pulling commode cushion, anti overturn device, the rear handle link the brakes, backrest can adjust the different angle, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, with safety belt.

66	52.5	44.5	45	20.5	43	118	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	26	92x26x95			



镀铬车架，高靠背，半躺，木质座厕板，附无孔PVC座垫，圆形灰色有盖深马桶，靠背可按需要调节成不同的角度，可拆卸高抬搁脚，8" 加厚前轮，24" 充气后轮，高强度塑料脚踏板，附安全带。

Chromed frame, high back, half lying, wood commode board, with imperforated PVC seat, rounded gray covered deep toilet, backrest can adjust the different angle, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, high strength plastic pedals, with safety belt.

64	52	39	45	25.5	44	104	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	26	92x29x96			



## 钢质手动轮椅系列 Steel Manual Wheelchair Series



**BWHM-1A905** 多功能型轮椅 Multifunctional wheelchair

镀铬车架，双交叉杆结构，高靠背，半躺，改进型可抽拉座厕垫，防后翻装置，后把手连刹车，靠背可按需要调节成不同的角度，可拆卸高抬搁脚，8"加厚前轮，24"充气后轮，带后刹车，附安全带，送餐桌板。

Chromed frame, double cross bar structure, high back, half lying, improved pulling commode cushion, anti overturn device, the rear handle link the brakes, backrest can adjust the different angle, detachable high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, with the rear brake, with safety belt and table board.

66	52.5	45.5	45	20.5	43	118	118
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	30	26	97x35x88			

## 铝合金轮椅系列 Aluminum Wheelchair Series



**BWHM-1A11(座厕垫附餐桌板)**

多功能型轮椅 Multifunctional wheelchair

铝合金阳极氧化车架，双交叉杆结构，固定扶手和搁脚，可抽拉式座厕垫，牛津布座靠垫，铝质小轮叉，8"实芯前轮，24"铝合金充气后轮，带后刹车，附安全带，送餐桌板。

Aluminum anodized frame, double cross bar structure, fixed armrest and footrests, pulling commode cushion, oxford cloth seat and backrest cushion, steel wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with the rear brake, with safety belt and table board.

64	42	41.5	41.5	21.5	43	104	87
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	13.5	92x29x88			

**BWHM-1A11** 普通轮椅 General wheelchair

## 铝合金轮椅系列 Aluminum Wheelchair Series



**BWHM-1A11** 普通轮椅 General wheelchair

铝合金阳极氧化车架，双交叉杆结构，固定扶手和搁脚，牛津布座靠垫，钢质小轮叉，8"实芯前轮，24"铝合金充气后轮，带后刹车，附安全带。

Aluminum anodized frame, double cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, steel wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with the rear brake, with safety belt.

64	42	41.5	41.5	21.5	43	104	87
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	13.5	92x29x88			



铝合金阳极氧化车架，双交叉杆结构，固定扶手和搁脚，牛津布座靠垫，钢质小轮叉，8"实芯前轮，24"铝合金充气后轮，带后刹车，附安全带。

Aluminum anodized frame, double cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, steel wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with the rear brake, with safety belt.

64	42	41.5	41.5	21.5	43	104	87
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	13.5	92x29x88			



## 铝合金轮椅系列 Aluminum Wheelchair Series

### BWHM-1A11 普通轮椅 General wheelchair



铝合金阳极氧化车架，双交叉杆结构，固定扶手和搁脚，牛津布座靠垫，钢质小轮叉，8" 实芯前轮，24" 铝合金充气后轮，带后刹车，附安全带。

Aluminum anodized frame, double cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, steel wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with the rear brake, with safety belt.

64	42	41.5	41.5	21.5	43	104	87
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	13.5	92x29x88			

### BWHM-1B33 便携式轮椅 Portable wheelchair



铝合金喷塑车架，双交叉杆结构，固定扶手和搁脚，牛津布座靠垫，铝合金小轮叉，6" 实芯前轮，14" 铝合金充气后轮，全车仅重10kg，附安全带。

Aluminum plastic-spray frame, double cross bar structure, fixed armrest and footrests, oxford cloth seat and backrest cushion, aluminum wheel fork, 6" solid front wheel, 14" aluminum pneumatic rear wheel, the whole wheelchair weighs only 10kg, with safety belt.

60	48	37	42	22.5	41	80	85
			重量 Weight(kg)	体积 Measurement(cm)			
30	15	26	12.5	83x24x87			

### BWHM-1B34 便携式轮椅 Portable wheelchair



### BWHM-1A11-1 普通轮椅 General wheelchair

铝合金阳极氧化车架，双交叉杆结构，固定扶手和搁脚，硬座，牛津布座靠垫，钢质小轮叉，8" 实芯前轮，24" 铝合金充气后轮，附安全带。

Aluminum anodized frame, double cross bar structure, fixed armrest and footrests, hard seat, oxford cloth seat and backrest cushion, steel wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with safety belt.

65	49	40	42.5	21.5	43	104	87
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	27	17	97x27x95			



铝合金阳极氧化车架，双交叉杆结构，可折背，固定扶手和搁脚，后把手连刹车，牛津布座靠垫，铝合金小轮叉，6" 实芯前轮，14" 铝合金充气后轮，带后刹车，全车仅重10kg，附安全带。

Aluminum anodized frame, double cross bar structure, folding back, fixed armrest and footrests, the rear handle link the brakes, oxford cloth seat and backrest cushion, aluminum wheel fork, 6" solid front wheel, 14" rear wheel, with the rear brake, the whole wheelchair weighs only 10kg, with safety belt.

60	48	37	42	22.5	41	80	85
			重量 Weight(kg)	体积 Measurement(cm)			
30	15	26	12.5	79x24x71			



## 铝合金轮椅系列 Aluminum Wheelchair Series

### BWHM-1B28 多功能型轮椅 Multifunctional wheelchair



铝合金阳极氧化车架，双交叉杆结构，翻转扶手，固定搁脚，牛津布座靠垫，后把手连刹车，铝合金小轮叉，8" 实芯前轮，24" 铝合金充气后轮，附安全带。

Aluminum anodized frame, double cross bar structure, turnover armrest, fixed footrests, oxford cloth seat and backrest cushion, the rear handle link the brakes, aluminum wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with safety belt.

70	46	41	44	20.5	46	103	87
			重量 Weight(kg)	体积 Measurement(cm)			
56	18	33	14	100x31x74			

### BWHM-1B28 多功能型轮椅 Multifunctional wheelchair



铝合金阳极氧化车架，双交叉杆结构，带后折背结构，翻转扶手，固定搁脚，牛津布座靠垫，后把手连刹车，铝合金小轮叉，8" 实芯前轮，24" 铝合金充气后轮，带后刹车，附安全带。

Aluminum anodized frame, double cross bar structure, with after folding back structure, turnover armrest, fixed footrests, oxford cloth seat and backrest cushion, the rear handle link the brakes, aluminum wheel fork, 8" solid front wheel, 24" aluminum pneumatic rear wheel, with the rear brake, with safety belt.

70	46	41	44	20.5	46	103	87
			重量 Weight(kg)	体积 Measurement(cm)			
56	18	33	14	100x31x74			

### BWHM-1A26 多功能型轮椅 Multifunctional wheelchair



### BWHM-1B28 多功能型轮椅 Multifunctional wheelchair

铝合金阳极氧化车架，双交叉杆结构，翻转扶手，固定搁脚，牛津布座靠垫，后把手连刹车，铝合金小轮叉，8" 实芯前轮，22" 铝合金免充气后轮，附安全带。

Aluminum anodized frame, double cross bar structure, turnover armrest, fixed footrests, oxford cloth seat and backrest cushion, the rear handle link the brakes, aluminum wheel fork, 8" solid front wheel, 22" aluminum free pneumatic rear wheel, with safety belt.

70	46	41	44	20.5	46	103	87
			重量 Weight(kg)	体积 Measurement(cm)			
56	18	33	14	100x31x74			



铝合金阳极氧化车架，双交叉管结构，翻转扶手，拆卸搁脚，牛津布座靠垫，8" 实芯前轮，24" 后轮，附安全带。

Aluminum anodized frame, double cross tube structure, turnover armrest, detachable footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" rear wheel, with safety belt.

64	51	41	42	22.5	45	103	89
			重量 Weight(kg)	体积 Measurement(cm)			
61	20	28	19	86x31x93			



## 铝合金轮椅系列 Aluminum Wheelchair Series

### BWHM-1A26 多功能型轮椅 Multifunctional wheelchair



铝合金阳极氧化车架，双交叉管结构，翻转扶手，拆卸搁脚，复式牛津布座靠垫，8" 实芯前轮，24" 后轮，带后刹车，附安全带。

Aluminum anodized frame, double cross tube structure, turnover armrest, detachable footrests, duplex oxford cloth seat and backrest cushion, 8" solid front wheel, 24" rear wheel, with the rear brake, with safety belt.

64	51	41	42	22.5	45	103	89

			重量 Weight(kg)	体积 Measurement(cm)
61	20	28	19	86x31x93

### BWHM-1B10 多功能型轮椅 Multifunctional wheelchair



铝合金阳极氧化车架，双交叉杆结构，高靠背，全躺，可拆卸高抬搁脚，8" 加厚前轮，24" 充气后轮，附安全带。

Aluminum anodized frame, double cross bar structure, high back, full lying, detachable high footrests, simple high footrests, 8" thicken front wheel, 24" pneumatic rear wheel, with safety belt.

66	52.5	44.5	45	20.5	43	118	118

			重量 Weight(kg)	体积 Measurement(cm)
61	20	30	26	92x26x96

### EA-6Q 旅行轮椅 Travel wheelchair

### BWHM-1A26 多功能型轮椅 Multifunctional wheelchair



铝合金阳极氧化车架，双交叉管结构，带后折背结构，翻转扶手，拆卸搁脚，牛津布座靠垫，8" 实芯前轮，24" 后轮，带后刹车，附安全带。

Aluminum anodized frame, double cross tube structure, with after folding back structure, turnover armrest, detachable footrests, oxford cloth seat and backrest cushion, 8" solid front wheel, 24" rear wheel, with the rear brake, with safety belt.

64	51	41	42	22.5	45	103	89

			重量 Weight(kg)	体积 Measurement(cm)
61	20	28	19	86x31x93



### 便携式多功能轮椅 Portable multifunctional wheelchair

可折叠便携式车架，实芯前后轮，自锁式后刹车，适合出行。

Portable and foldable frame, solid front and rear wheel, self-locking rear brake, suitable for travel.

51	45	40	41	25	38	41	89

			重量 Weight(kg)	体积 Measurement(cm)
15	12.5	21	6.9	46x34x66



## 电动轮椅系列 Electric Wheelchair Series



### BWHE-1A26D

老年电动代步车 The old electric vehicle controller

铝合金车架，高级，节能电池；万向，省力驱动；8" 发泡前轮，24" 后轮。

Aluminum frame, advanced, energy saving battery; universal, labor saving drive; 8" foam front wheel, 24" rear wheel



### BWHE803

老年电动代步车 The old electric vehicle controller

灵活、轻便、舒适，加厚柔软座靠垫，橡胶发泡减震前后轮。

Flexible, lightweight, comfortable, thickened soft seat and backrest cushion, foam rubber shock absorber front and rear wheel.

## 三轮车系列 3-wheel Rollator Series



### BWBM-346A-1 手摇式三轮车 Hand 3-wheel wheelchair

### BWBM-346A 推拉式三轮车 Push-pull 3-Wheel rollator

利用曲柄连杆运动原理，省掉链条、齿轮，单手轻轻推拉，即可实现前进、后退及转弯，非常适用于下肢瘫痪、行动不便者使用。

Exploiting movement principle of crank connecting link, chains and gears are no longer necessary. Push and pull with one hand slightly, to go forward and back and make a turn. It's very suitable for person with paraplegia of the lower limbs or limited physical mobility.

64	46	41	46	24	48	90	152
			重量 Weight(kg)	体积 Measurement(cm)			
90	38	38	31	78x67x64			



手摇三轮车轮椅可供下肢行动不便者乘坐，双手摇动摇把驱动前轮前进或制动。

Hand 3-wheel wheelchair available for lower limb disabled person ride, hands shaking handle driving the front wheels forward or brake.

64	46	41	46	24	48	90	152
			重量 Weight(kg)	体积 Measurement(cm)			
90	38	38	31	78x67x64			



## 手/拐/腋杖系列 Walking Stick / Crutch / Armpit Crutch Series



BWHM601

BWHM103

单脚拐，铝合金表面阳极处理，颜色可根据需要制作。手柄与脚垫的颜色相吻合，高度可调节，每调节一格高度增加1"。

Stick, constructed of anodized aluminum finish, color by order. Matching color hands and tip, height adjusts in 1" increments.

型号 Model NO.	高(cm) Height	净重(kg) N.W
BWHM601	89-120	0.5
BWHM103	68-90	0.3
BWHM801	84	1

### BWHM801 多功能手杖 Multifunctional crutch

带凳三脚拐，铝合金表面阳极处理，在行走累时可展开当座凳使用，收展自如，方便实用。

Three feet cane with stool, constructed of anodized aluminum finish, a comfortable seat when open, tired when walking can be expanded to use when stool, collapse and expand very easy, convenient and practical.



### BWHM161 腋拐 Crutch

铝制腋下拐，产品为铝材质表面阳极处理，非常轻便，泡棉护肩，泡棉握把为高密度的橡胶发泡而成，塑胶拐柄，握把为高强度的塑胶制成，是很好的护健产品。按不同人身高分为大、中、小型。

Aluminum crutches, constructed of anodized aluminum finish, very lightweight, completely assembled with crutch foam cushions, handgrips and reinforced metal rubber tips ready to fit, it is a very good product for health care. According to different height is divided into large, medium, small.



S

M

L

型号 Model NO.	高(cm) Height	净重(kg) N.W
BWHM161(L)	129-149	0.9
BWHM161(M)	116-136	0.85
BWHM161(S)	103-123	0.8



BWHM402

BWHM404

型号 Model NO.	高(cm) Height	净重(kg) N.W
BWHM402	71-94	0.9
BWHM404	69-91	0.9

### BWHM402 四脚拐 4-foot walking stick

大四脚拐，主杆为铝材质表面阳极处理，底座为铁板、铁管焊接而成；四个分支给一种平衡感；可调高度：29" ~38"。

The main bar constructed of anodized aluminum, the base for iron, iron pipes and welded; Offset handle provide balance for steady handing and feeling; Adjustable height 29" ~38" .

### BWHM404

四脚拐，主杆为铝材质表面阳极处理，底座为铁板、铁管焊接而成；四个分支给一种平衡感；可调高度：29" ~38"。

The main bar constructed of anodized aluminum, the base for iron, iron pipes and welded; Offset handle provide balance for steady handing and feeling; Adjustable height 29" ~38" .



## 助行器系列 Walking Aid Series

## 座厕椅系列 Commode Chair Series



尺寸/Size(cm): 56(L)x56(W)x(74.5-92)(H)



尺寸/Size(cm): 57(L)x60(W)x(74.5-92)(H)



### BWHM-1B16 座厕椅 Commode chair

加强型座厕椅，托架、靠背可以折叠，携带方便，不需任何工具即可组装，高度可调。

Enhanced commode chair, the bracket, backrest and armrest are all foldable, so the freight will be saved, no extra tools are needed for assembling, adjustable in height.

### BWHM-407 助行器 Walking aid

适用于老年人，腿、脚行动不方便者出外散步使用；产品主体为铝材质，比较轻便，并可折叠，不占地方。

Applicable to the elderly or leg, foot action inconvenient to go out for a walk; The main body of the product is lead material, more portable, and can be folded, do not take place.

### BWHM-408 助行器 Walking aid

此助行器前面一对1/2"小轮，使用起来更方便、更安全。

In front of the walking aid with one pair 1/2" wheels, to use of more convenient, more secure.



尺寸/Size(cm): 59(L)x50(W)x(83-91)(H)



尺寸/Size(cm): 73(L)x68(W)x(83-92)(H)



### BWHM-1B23 座厕椅 Commode chair

加强型座厕椅，托架、靠背可以折叠，携带方便，不需任何工具即可组装，高度可调。

Enhanced commode chair, the bracket, backrest and armrest are all foldable, so the freight will be saved, no extra tools are needed for assembling, adjustable in height.



### BWHM-1B7 座厕椅 Commode chair

固定扶手座厕椅，宜室内使用。

Fixed armrest commode chair, suitable to use indoor.

### BWHM-1B21 助行器 Walking aid

可折叠，适用于超市购物，喷塑车架，带座垫，手柄高度可调节，带线刹，附网篮和托盘，8" 聚氨酯前后轮。

It can be folded, suit for shopping, plastic-spray frame, with seat, adjustable handle height, with cable brake, with baskets and trays, 8" PU front and rear wheels.

### BWHM-405 助行器 Walking aid

折叠式两轮助行器，钢管镀铬，手柄高度可调节，带座垫。

Foldable two wheels walking aid, chromed steel frame, adjustable handle height, with seat.

型号 Model NO.										重量 Weight(kg)	体积 Measurement(cm)
BWHM-1B7	46	50	42	20	20	48	50	91	48	10	58x21x92