

## General Description:

WANROOEMED **Patient Lift** is used for reducing the possibility of caregiver back injury and to ensure dignity in patient handling. These high-quality patient lifts are comfortable, reliable and integral tools for staff and caregivers. Using the patient lift for patient handling can help reduce staff injuries, save you money and reinforce safety as a priority in your health care facility

## DG201 Battery-Powered Patient Lift



## Features

- Leg base open easily by foot pedal
- Height range allows lifting from floor and high surfaces.
- manual emergency stop for safety
- Long, padded handles offer a plethora of grip choices.
- Patient can turn 360° at any height during the lifting, to access a most comfort position.
- The hoister handle design accords to dynamics, easy to push forward, backward, or turn.
- Lower base height maximizes stability and fits under most

beds

- Warning beep for low battery charge and weight capacity levels

## Specifications

|                               |  |
|-------------------------------|--|
| <b>Weight</b>                 | 43kgs                                  |
| <b>Weight Capacity</b>        | 150kgs standard <b>200kgs optional</b> |
| <b>Overall Length</b>         | 1100mm                                 |
| <b>Lifting range</b>          | 740-1775mm                             |
| <b>Legs open width</b>        | 800mm                                  |
| <b>Legs closed width</b>      | 585mm                                  |
| <b>Overall height of legs</b> |  |
| <b>Ground clearance</b>       | 120mm                                  |
| <b>Front twin casters</b>     | 2"                                     |
| <b>Rear break casters</b>     | 4" with brakes                         |
| <b>Power</b>                  | 220V,50HZ standard                     |
| <b>Battery</b>                | 24V Rechargeable sealed                |
| <b>Charger</b>                | Output: 24V DC                         |

no order

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com)     TEL: 0086-512-58810860     FAX: 086-512-58810663

|                 |  |
|-----------------|--|
| <b>Warranty</b> | Limited One (1) Year Warranty  |
| <b>Package</b>  | 1160*660*250mm, G.W.: 48kgs,N.W.: 43kg,156pcs/20GP, 306pcs/40GP, 346pcs/40HQ |

## DG202 Battery-Powered Stand-up Lift

### DG202 Battery-Powered Stand-up Lift

WANROOEMED DG202 Battery-Powered Stand-up Lift is one ideal transport device to provide quick, safe assistance to those who have trouble standing. With multi-safety design, DG202 Battery-Powered Stand-up Lift ensures caregiver and resident security

#### Features

- Heavy-gauge steel construction
- Leg open function to mostly close the wheel chair or medical bed.
- Height range allows lifting from floor and high surfaces.
- manual emergency stop for added safety
- Adjustable knee pad unit for comfort at any height.
- 24V DC motor allows quiet, smooth operation
- Warning beep for low battery charge and weight capacity levels
- Lower base height maximizes stability and fits under most beds



#### Specifications

|                        |                 |                        |
|------------------------|-----------------|------------------------|
| <b>Weight Capacity</b> | 150kgs standard | <b>200kgs optional</b> |
|------------------------|-----------------|------------------------|

|                       |        |
|-----------------------|--------|
| <b>Overall Length</b> | 1200mm |
|-----------------------|--------|

|                      |            |
|----------------------|------------|
| <b>Lifting range</b> | 760-1820mm |
|----------------------|------------|

|                        |       |
|------------------------|-------|
| <b>Legs open width</b> | 810mm |
|------------------------|-------|

|                          |       |
|--------------------------|-------|
| <b>Legs closed width</b> | 640mm |
|--------------------------|-------|

|                               |  |
|-------------------------------|--|
| <b>Overall height of legs</b> |  |
|-------------------------------|--|

|                         |       |
|-------------------------|-------|
| <b>Ground clearance</b> | 120mm |
|-------------------------|-------|

|                           |    |
|---------------------------|----|
| <b>Front twin casters</b> | 2" |
|---------------------------|----|

|                           |                |
|---------------------------|----------------|
| <b>Rear break casters</b> | 4" with brakes |
|---------------------------|----------------|

|              |                    |
|--------------|--------------------|
| <b>Power</b> | 220V,50HZ standard |
|--------------|--------------------|

|                |                         |
|----------------|-------------------------|
| <b>Battery</b> | 24V Rechargeable sealed |
|----------------|-------------------------|

no order

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 FAX: 086-512-58810663

**Charger** Output: 24V DC

**Warranty** Limited One (1) Year Warranty

**Package** 1160\*660\*300mm,G.W.: 53kgs,N.W.: 48kgs 150pcs/20GP, 300pcs/40GP,

## DG203 Battery-Powered Assist Lift

### DG203 Battery-Powered Assist Lift

WANROOEMED Battery-Powered Assist Lift is one ideal transport device to provide quick, safe assistance to those who have trouble standing. With multi-safety design, DG203 Battery-Powered Stand-up Lift ensures caregiver and resident security

#### Features

- Heavy-gauge steel construction
- Leg open function to mostly close the wheel chair or medical bed.
- Height range allows lifting from floor and high surfaces.
- manual emergency stop for added safety
- Adjustable knee pad unit for comfort at any height.
- 24V DC motor allows quiet, smooth operation
- Warning beep for low battery charge and weight capacity levels
- Lower base height maximizes stability and fits under most beds



#### Specifications

**Weight Capacity** 150kgs standard **200kgs optional**

**Overall Length** 1100mm

**Lifting range** 390-1480mm/680-1280mm

**Legs open width** 840mm

**Legs closed width** 650mm

**Overall height of legs**

no order

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 FAX: 086-512-58810663


|                           |  |
|---------------------------|--|
| <b>Ground clearance</b>   | 120mm  |
| <b>Front twin casters</b> | 2"   |
| <b>Rear break casters</b> | 4" with brakes   |
| <b>motor</b>              | 220V,50HZ standard,8000N   |
| <b>Battery</b>            | 24V Rechargeable sealed  |
| <b>Charger</b>            | Output: 24V DC   |
| <hr/>                     |  |
| <b>Warranty</b>           | Limited One (1) Year Warranty  |
| <b>Package</b>            | 1160*660*300mm,G.W.: 60kgs,N.W.: 55kgs 150pcs/20GP, 300pcs/40GP, 340pcs/40HQ |

## DG100 Hydraulic Patient Hoist

### DG100 Hydraulic Patient Hoist



DG100 Hydraulic Patient Hoist is suitable for stay in bed for a long term, Patients with inconvenient action or can not move after operation. They use it accompanied by nursing staff, Hospital, nursing center or family can be used. The main frame is made of high strength steel, the powder of surface plastic spraying treatment, The front wheel of the castor is 3 inch universal wheel, the back wheel is 4 inch brake wheel, foot stool can adjust leg width, Sling suspension is safe and reliable, The hydraulic lifting device, the scope is broader.



Lifting handle



Adjustable foot pedal



Release handel

*no order*

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 FAX: 086-512-58810663

## Specifications

|                           |   |
|---------------------------|---|
| <b>Weight Capacity</b>    | 150kgs standard   |
| <b>Overall Length</b>     | 1100mm  |
| <b>Pull handle height</b> | 1130mm  |
| <b>Lifting range</b>      | 740-1775mm  |
| <b>Legs open width</b>    | 800mm   |
| <b>Ground clearance</b>   | 120mm   |
| <b>Front twin casters</b> | 3"  |
| <b>Rear break casters</b> | 4" with brakes  |
| <b>Warranty</b>           | Limited One (1) Year Warranty   |
| <b>Package</b>            | 1160*660*250mm, 144pcs/20GP, 324pcs/40GP, 360pcs/40HQ, gross weight: 41kgs, net weight: 38kgs |

## Patient Lift Slings

### DGA101 Economy All-purpose Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

#### features

- green nylon knitted fabric
- easy to install, adjust
- excellent back support and protect
- gray-blue
- machine washable
- loading 150KG



### DGA406NE Standard Grid Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

#### features

no order

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 FAX: 086-512-58810663

- green nylon knitted fabric
- easy to install, adjust
- excellent back support and protect
- gray-blue
- machine washable
- loading 150KG
- It has large, medium and small size

## DGA407NE Standard Mesh Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

### features

- green nylon knitted fabric or green nylon mesh cloth
- easy to install, adjust
- excellent back and head support and protect
- fit to buttocks area
- leg pad
- gray-blue
- machine washable
- loading 150KG
- It has large, medium and small size



## Patient Lift Slings

## DGA400SF Standing Type Patient Lift Sling

- This product is suitable for patient lift of family, hospital, sanatorium etc.
- features
- green nylon knitted fabric
- easy to connect and adjust
- protect back support, it has fixed belt
- suitable for standing type or multi-functional type patient lift
- gray-blue
- 150KG; loading 150KG
- machine washable
- It has large, medium and small size



## DGA406N Multi-functional Type Patient Lift Sling

- This product is suitable for patient lift of family, hospital, sanatorium etc.
- features
- green nylon knitted fabric or green mesh fabric
- easy to fit and adjust
- not fit to buttocks and head area, excellent back protection
- leg pad



no order

All specifications, pictures are subject to change with  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 F



- it has fixed belt
- gray-blue
- loading 150KG
- machine washable
- It has large, medium and small size

### DGA409N Hole Type Patient Lift Sling

- This product is suitable for patient lift of family, hospital, sanatorium etc.
- features
- green nylon knitted fabric or green nylon mesh cloth
- easy to fit and adjust
- excellent back and head support and protect
- with toilet hole
- gray-blue
- loading 150KG
- machine washable
- It has large, medium and small size



## Patient Lift Slings

### DGA408N Whole Body Protection Type Patient Lift Sling

This product is suitable for patient lift of family, hospital, sanatorium etc.

features

- green nylon knitted fabric or green mesh fabric
- easy to fit and adjust
- whole body protection include head support and legs can be divided
- suitable for the patient who needs whole body nursing
- with toilet seat hole
- leg pad
- gray-blue
- machine washable
- loading: 150KG
- It has large, medium and small size



**NOTICE: N type made of nylon knitting material, NE type made of nylon mesh material. All available according to requirement.**

no order

All specifications, pictures are subject to change without earlier notification.  
[www.wanrooemed.com](http://www.wanrooemed.com) TEL: 0086-512-58810860 FAX: 086-512-58810663