



HBHKN-9010

INFANT RADIANT WARMER



Stock code [300314]



JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.



Brief Introduction

- Three control Modes: Pre-warm Mode, Manual Mode and Baby Mode control by micro-computer;
- Set Temperature and Baby Temperature can be displayed separately;
- Horizontal angle of warmer module and inclination of the bassinet can be adjusted;
- The panels surrounded bassinet can be turned outward or disassembled;
- X-ray cassette under the bassinet;
- APGAR timer;
- Multiple failure alarm indications;
- RS232 connector for data output;
- Disposable skin temperature sensor and SpO₂ are optional.

Standard Configuration

Main body (including the warmer module, bassinet, main column, mobile stand and temperature controller), I.V. Pole and tray.

Specification

Power requirement: AC220V-230V/50Hz or AC110-120V/50-60Hz, 600VA

Control Mode: Pre-warm Mode, Manual Mode and Baby Mode are controlled by micro-computer

Baby Mode temperature control range: 34.5-37.5°C

Temperature sensor display range: 5°C-65°C

Accuracy of skin temperature sensor: $\pm 0.2^{\circ}\text{C}$

Temperature uniformity of mattress: $\leq 2^{\circ}\text{C}$

Angle of heater head: 0° -90° any angel for two ways

Inclination of the bassinet: three tilt position

APGAR timer: audible and visible alarm when the device runs to 50" ~ 1' , 4' 50" ~ 5' , 9' 50" ~ 10'

Failure alarm: over temperature alarm, deviation alarm, sensor failure alarm, power failure alarm, setting alarm, checking alarm and so on.

Environmental Requirement

Operating range: 18°C~30°C

Ambient air velocity rate: <0.3m/s

Product Packaging

Package: The whole unit is packed into one carton

Dimension: 1220mm x 670mm x 1100mm

Gross Weight: 100kg



JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19, East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663 Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanrooemed.com Website: www.rooemed.com www.wanrooemed.com