

USC60

LAPTOP COLOR DOPPLER
ULTRASOUND SCANNER (BASIC
4D MODE) /

B, 2B, 4B, BM, B/C, B/D B/C/D, B/CFM/D, Color, Dual Color PW, Duplex/ Triplex, CFM, PD, Directional PD real time 4D



Benefits

Powerful, professional measurement tools allow you to give an accurate and quick diagnosis. The measurement packages of the ultrasound are comprehensively applicable to cardiology,OB/GYN,radiology,anesthesiology,musculoskeletal department etc. The high-level ultrasound diagnosis systems provide excellent solutions for both adults and pediatrics. You can make clinical decisions with greater confidence and expand your practice in new directions.

Features

15-inch high-definition medical LED monitor

Compact and streamlined design

PW spectrum with measurement

CFM, PDI, DPDI, HPRF, THI and etc.

Dual/Triplex synchromous technology


Buit-in lithium battery (more than 3 hours standby time)

Support multiple data transmission interface

4D imaging

Optional trolley and suitcase

3 USB Port, HDMI, Dicom 3.0 supported



High-density and broadband frequency transducers

