

USMSU3

Full-digital Mechanical Sector Ultrasound Scanner



Pregnancy test



Backfat measurement



Compact and
lightweight



Large capacity
battery




USMSU3

Full-digital Mechanical Sector Ultrasound Scanner



- Latest cortex-m4 embedded control system + mini-ultrasound technology hardware system. A scan and B scan
- 3.5" high resolution color LCD, smart, light and portable, suitable for outdoor use
- Measure back fat thickness, loin muscle thickness and lean percentage
- 3 kinds scale of image display
- With U disk, image and measured data can be stored in Excel table, easy for data management processing
- Unique character hidden function, automatically hidden in the setting time
- System automatically freeze and protect probe with time setting
- Customized languages: English, French, German, Spanish, Portuguese, Russian, etc.
- Lithium battery: Large capacity (2600mAh), working time $\geq$ 300 minutes
- Net weight: 620g



万容机械有限公司 WANROOE MACHINERY CO., LTD.
 江苏省张家港市人民东路19号国泰新天地广场1311,1313室
 Room1311,1313, Guotai New World Plaza, No19, East Renmin Road,Zhangjiagang,
 Jiangsu, China.
 Tel:+86 512 58810663 Fax:+ 86 512 58810663
 Email:doris@wanrooemed.com