

support you care you



WANROOEMED

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY

AOT700A

AOT700A



support you care you

AOT700A

WELL SOLD FOR 40 YEARS WORLDWIDE
MEET THE BASIC DEMANDS OF GENERAL SURGERY



Introduction

700A is a multi-functional operating table that meets various requirements of operations. long horizontal sliding, suitable for C arm and X-ray. The "T-shape" base design makes enough free space for the doctors. High automation, low noise, high reliability. Key parts imported from abroad, can be regarded as an ideal electric table.

Main Technical Data

Parameters may have reasonable error.

Item	AOT700A
Length and width for tabletop	1970/500mm
Elevation(down/up)	690/1100mm
Trendelenburg/anti-trendelenburg	30°/30°
Lateral tilt(left/right)	18°/18°
Head plate(up/down)	60°/90°
Leg Plate(up/down/outward)	15°/90°/90°
Back plate(up/down)	90°/45°
Horizontal sliding	350mm
Power	220V , 50Hz, 1.0kW

江苏如翼医疗科技有限公司

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19, East Renmin Road, Zhangjiagang, Jiangsu, China.

Tel: +86 512 58810663 Whatsapp/Mobile: + 86-18962287398; 15962399164

Email: info@wanrooemed.com Http://wanrooemed.com