


F-PS60L Plasma Sterilizer



Features:

- 1.7 inch LCD(liquid crystal display) and Android operating system;
- 2.Support wired LAN,WLAN and 3G;
- 3.2 USB connectors,which can export data record;
- 4.Text message reminder;
- 5.Video broadcast;
- 6.Cloud service,backup and recovery record data;
- 7.Real time curve:the curve shows the sterilization status for easy monitoring;
- 8.Online infection control knowledge,user manual and maintenance instructions are available;
- 9.update online medical industry news;
- 10.Update application online and offline;
- 11.Export record data as excel file;
- 12.Leybold vacuum pump from Germany.



to order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

| | |
|--------------------|--|
| Type | F-PS60 |
| Overall Dimension | H1725*W750*D900(mm) |
| Door Type | Single automatic vertical sliding door |
| Capacity | 60L(cylinder chamber) |
| Chamber Size | Φ450*690(mm) |
| Sterilization Mode | Quick/standard/enhanced cycle |
| Sterilization Temp | 50°C--60°C |
| Power Rating | 220V/50Hz,2.7Kw |
| Control system | Android system,intelligent management |
| Cassette | 12 capsules for each cassette which can be used for 6 cycles |

to order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663