

ISO9001
ISO13485



A-C102系列多功能产床

A-C102 MULTIFUNCTION OBSTETRIC TABLES

20年专注妇产科设备研发生产



万容机械有限公司
WANROOE MACHINERY CO., LTD.

张家港市金港大道288号辰龙大厦

Chenlong Building, No 288, Jingang Road,
Zhangjiagang, Suzhou, Jiangsu, China

Tel: +86 512 58810860

Fax: +86 512 58810663

Email: info@wanrooemed.com



[Http://wanrooemed.com](http://wanrooemed.com)



中国·江苏
JIANGSU CHINA

张家港市万容机械有限公司
WANROOE MACHINERY CO., LTD.

A-C102 多功能产床

A-C102 MULTIFUNCTION OBSTETRIC TABLE



床面长度及宽度	长度1920mm,宽度610mm
床面最低最高高度	790mm ~ 1040mm
背板折转角度	-5° ~ 75°
座板上折角度	0° ~ 15°
辅助台长度及宽度	长度700mm,宽度610mm
辅助台外摆角度	0° ~ 90°
搁手板外摆角度	≥90° (可向上折起作为护栏使用)
污物盆尺寸	468 × 338 × 125mm (可带排污软管)
电源	AC 220V 50HZ

Table top dimension: length 1920mm, width 610mm
Min. & Max. table height: 790mm ~ 1040mm
Back section turning: -5° ~ 75°
Seat section turn up: 0° ~ 15°
Assist platform dimension: length 700mm, width 610mm
Assist platform swing out: 0° ~ 90°
Arm rest swing out: ≥90° (can be at upright position as guard rail)
Filth basin dimension: 468 × 338 × 125mm (can be equipped with soft pipe)
Power: AC 220V 50HZ

附 件	Standard accessories
防水垫 1件	Waterproof cushion 1 set
腰部扶手 1付	Waist grab handle 1 pair
搁腿架 1付	Leg holder 1 pair
污物盆 1只	Filth basin 1 piece
搁手板 1付	Arm rest 1 pair
脚踏开关 1只	Foot switch 1 piece
承合器 1付	Clamp 1 pair

备注：床面为无缝皮革,模具成形
Note: Mattress made from mould, with no seam on the surface.



A-C102A 多功能产床

A-C102A MULTIFUNCTION OBSTETRIC TABLE



床面长度及宽度	长度1900mm,宽度610mm	Table top dimension: length 1900mm, width 610mm
床面最低最高高度	750mm ~ 1000mm	Min. & Max. table height: 750mm ~ 1000mm
床面前后倾角度	前倾 $\geq 15^{\circ}$, 后倾 $\geq 15^{\circ}$	Trendelenburg & reverse trendelenburg: reverse trend $\geq 15^{\circ}$, trend $\geq 15^{\circ}$
背板折转角度	$0^{\circ} \sim 75^{\circ}$	Back section turning: $0^{\circ} \sim 75^{\circ}$
辅助台长度及宽度	长度550mm,宽度550mm	Assist platform dimension: length 550mm, width 550mm
搁手板外摆角度	$\geq 90^{\circ}$	Arm rest swing out: $\geq 90^{\circ}$
污物盆尺寸	468 $\times$ 338 $\times$ 125mm (可带排污软管)	Filth basin dimension: 468 $\times$ 338 $\times$ 125mm (can be equipped with soft pipe)
电源	AC 220V 50HZ	Power: AC 220V 50HZ

附件	Standard accessories	
防水垫 1件	Waterproof cushion	1 set
腰部扶手 1付	Waist grab handle	1 pair
搁腿架 1付	Leg holder	1 pair
污物盆 1只	Filth basin	1 piece
搁手板 1付	Arm rest	1 pair
脚踏开关 1只	Foot switch	1 piece
承合器 1付	Clamp	1 pair

备注: 床面为无缝皮革,模具成形
Note: Mattress made from mould, with no seam on the surface.

A-C102B 多功能产床

A-C102B MULTIFUNCTION OBSTETRIC TABLE



床面长度及宽度	长度1900mm,宽度610mm	Table top dimension: length 1900mm, width 610mm
床面最低最高高度	750mm ~ 1000mm	Min. & Max. table height: 750mm ~ 1000mm
床面前后倾角度	前倾 $\geq 15^{\circ}$, 后倾 $\geq 15^{\circ}$	Trendelenburg & reverse trendelenburg: reverse trend $\geq 15^{\circ}$, trend $\geq 15^{\circ}$
背板折转角度	$0^{\circ} \sim 75^{\circ}$	Back section turning: $0^{\circ} \sim 75^{\circ}$
辅助台长度及宽度	长度550mm,宽度550mm	Assist platform dimension: length 550mm, width 550mm
搁手板外摆角度	$\geq 90^{\circ}$	Arm rest swing out: $\geq 90^{\circ}$
脚踏外摆角度	$0^{\circ} \sim 45^{\circ}$	Footrest swing out: $0^{\circ} \sim 45^{\circ}$
脚踏上折角度	$0^{\circ} \sim 45^{\circ}$	Footrest turn up: $0^{\circ} \sim 45^{\circ}$
污物盆尺寸	468 $\times$ 338 $\times$ 125mm (可带排污软管)	Filth basin dimension: 468 $\times$ 338 $\times$ 125mm (can be equipped with soft pipe)
电源	AC 220V 50HZ	Power: AC 220V 50HZ

附件	Standard accessories	
防水垫 1件	Waterproof cushion	1 set
腰部扶手 1付	Waist grab handle	1 pair
搁腿架 1付	Leg holder	1 pair
污物盆 1只	Filth basin	1 piece
搁手板 1付	Arm rest	1 pair
脚踏开关 1只	Foot switch	1 piece
承合器 1付	Clamp	1 pair

备注: 床面为无缝皮革,模具成形
Note: Mattress made from mould, with no seam on the surface.



A-C102C 多功能产床

A-C102C MULTIFUNCTION OBSTETRIC TABLE



床面长度及宽度	长度1900mm,宽度610mm
床面最低最高高度	660mm~910mm
床面前后倾角度	前倾 $\geq 10^{\circ}$,后倾 $\geq 15^{\circ}$
背板折转角度	$-10^{\circ} \sim 75^{\circ}$
辅助台长度及宽度	长度550mm,宽度550mm
搁手板外摆角度	$\geq 90^{\circ}$
污物盆尺寸	468 $\times$ 338 $\times$ 125mm (可带排污软管)
电源	AC 220V 50HZ

Table top dimension: length 1900mm, width 610mm
 Min. & Max. table height: 660mm ~ 910mm
 Trendelenburg & reverse trendelenburg: reverse trend $\geq 10^{\circ}$, trend $\geq 15^{\circ}$
 Back section turning: $-10^{\circ} \sim 75^{\circ}$
 Assist platform dimension: length 550mm, width 550mm
 Arm rest swing out: $\geq 90^{\circ}$
 Filth basin dimension: 468 $\times$ 338 $\times$ 125mm (can be equipped with soft pipe)
 Power: AC 220V 50HZ

附 件	Standard accessories
防水垫 1件	Waterproof cushion 1 set
腰部扶手 1付	Waist grab handle 1 pair
搁腿架 1付	Leg holder 1 pair
污物盆 1只	Filth basin 1 piece
搁手板 1付	Arm rest 1 pair
脚踏开关 1只	Foot switch 1 piece
承合器 1付	Clamp 1 pair

备注: 床面为无缝皮革,模具成形
 Note: Mattress made from mould, with no seam on the surface.



A-C102D-01 多功能产床 (脚踩液压)

A-C102D-01 MULTIFUNCTION OBSTETRIC TABLE (HYDRAULIC)

床面长度及宽度	长度1880mm,宽度610mm
床面最低最高高度	670mm~920mm
床面前后倾角度	前倾 $\geq 10^\circ$,后倾 $\geq 20^\circ$
背板折转角度	$-6^\circ \sim 70^\circ$
辅助台长度及宽度	长度540mm,宽度610mm
搁手板外摆角度	$\geq 90^\circ$
ABS污物盆容积	14升

Table top dimension: length 1880mm, width 610mm
Min. & Max. table height: 670mm ~ 920mm
Trendelenburg & reverse trendelenburg: reverse trend $\geq 10^\circ$, trend $\geq 20^\circ$
Back section turning: $-6^\circ \sim 70^\circ$
Assist platform dimension: length 540mm, width 610mm
Arm rest swing out: $\geq 90^\circ$
ABS filth basin capacity: 14L

附件	Standard accessories
防水垫 1件	Waterproof cushion 1 set
腰部扶手 1付	Waist grab handle 1 pair
搁腿架 1付	Leg holder 1 pair
污物盆 1只	Filth basin 1 piece
搁手板 1付	Arm rest 1 pair
承合器 1付	Clamp 1 pair

备注: 床面为无缝皮革,模具成形
Note: Mattress made from mould, with no seam on the surface.



标准附件

Standard accessories



腰部扶手
Waist grab handle



承合器
Clamp



防水垫
Waterproof cushion



搁手架
Arm rest

选配件 Optional Accessories



搁腿架
Leg holder



带软垫搁腿架
Leg holder with soft pad



不锈钢污物盆
Stainless steel filth basin



工程塑料污物盆
Engineering plastic filth basin



手控器
Hand control



脚控器
Foot switch



带脚蹬搁腿架
Leg holder with footrest



脚蹬
Footrest

附加件 Additional Accessories



移动式污物桶
Kick bucket



头部拉手
Head section handle



可移动辅助台
Moveable assist platform



可升降辅助台
Height adjustable assist platform

A-C102D-02 多功能产床 (脚踩液压)

A-C102D-02 MULTIFUNCTION OBSTETRIC TABLE (HYDRAULIC)

床面长度及宽度	长度1895mm,宽度610mm
床面最低最高高度	630mm~880mm
床面前后倾角度	前倾 $\geq 10^\circ$,后倾 $\geq 15^\circ$
背板折转角度	$-10^\circ \sim 75^\circ$
辅助台长度及宽度	长度550mm,宽度610mm
搁手板外摆角度	$\geq 90^\circ$
污物盆尺寸	468×338×125mm (可带排污软管)

Table top dimension: length 1895mm, width 610mm
Min. & Max. table height: 630mm ~ 880mm
Trendelenburg & reverse trendelenburg: reverse trend $\geq 10^\circ$, trend $\geq 15^\circ$
Back section turning: $-10^\circ \sim 75^\circ$
Assist platform dimension: length 550mm, width 610mm
Arm rest swing out: $\geq 90^\circ$
Filth basin dimension: 468×338×125mm (can be equipped with soft pipe)

附件	Standard accessories
防水垫 1件	Waterproof cushion 1 set
腰部扶手 1付	Waist grab handle 1 pair
搁腿架 1付	Leg holder 1 pair
污物盆 1只	Filth basin 1 piece
搁手板 1付	Arm rest 1 pair
承合器 1付	Clamp 1 pair

备注: 床面为无缝皮革,模具成形
Note: Mattress made from mould, with no seam on the surface.

