

ITEM		TECHNICAL REQUIREMENT
Size	Length	1930mm
	Width	645mm
Range Of Motion	Back lift	0~70°
	High-low	510-850mm
Main Material	Bed Board, Frame	Bed board-Integrated molding of PP resin. Frame-Aluminum alloy pressing and high quality steel.
	Mattress	Surface material-Polyester Inner core material-polyurethane foam
Diameter of Caster	Double sided casters with central lock	150mm
Safe Working Load		170kg



Standard Function

- ✦ Bed Body
- ✦ Movable Mattress
- ✦ Side Rails
- ✦ Double Sided Caster
- ✦ Central Lock
- ✦ Fifth Wheel
- ✦ Base Cover
- ✦ Oxygen Bottle Holder
- ✦ Hand Crank
- ✦ IV Pole



BD111-3

Applicable Department: Emergency room /
Gastroscope room / Operation room
Multifunctional Stretcher Cart

CARE AND LOVE ALWAYS WITH YOU



Rotating Side Rails

- ✦ Side rails can be fixed at a horizontal position for drip and puncture . Loading capacity 10kg.
- ✦ Concave design can prevent catheter slide.



Aluminum Alloy Pressing

- ✦ Integrated molding, more strength , more style.



Back Lift

- ✦ Using gas spring to Stepless lift.



Functional Presentation Of Bed Body

- ✦ Back up and down 0-70 degrees , high-low 510-850mm.



Oxygen Bottle Holder

- ✦ When standby, put it under the backboard, easy and safe.
- ✦ Oxygen capacity 500 liters (Diameter 105-115mm).



Movable Mattress

- ✦ Using high tech waterproof fabric and electro-static prevention, Washable and easy clean, design of three sections, only one person can transfer the patient.



Silent Casters With Central Lock

- ✦ 150mm diameter resin casters with lock pedal on four Corners, easy for nurse to operate.



Fifth Wheel

- ✦ The conversion of stretcher cart is easily realized between "straight" and "free" by operating the lever.
- ✦ Easier to control the direction with "straight".



Base Cover

- ✦ Base cover has two sections of different size and depth, multiple leaking holes, easy to use and clean, loading capacity 10kg.



Double Locks Of Side Rails

- ✦ Double lock on the foot side, prevent wrong operation, more safe.

