

3-Function Manual up/down emergency stretcher BD207-1



Back Rest Tilting $\geq 70^\circ$

Footrest Tilting

Hi-Low

Features:

- 1)The emergency bed is used for carrying patients and the wounded people for the hospitals
- 2)Material: high quality stainless steel and polishing
- 4)which could make the operation more relaxing
- 5)Four 5" silent castors.two with break
- 6)1 pair of stainless steel side rails
- 7)3 sets of cranks

to order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663

8)1 piece of shelf

9)1 piece of Oxygen cart holder

10)2 Urine hooks

11)ABS foldable crank handles

12)Switch of the side rails

Dimentions:

Size:1950*700*450~680mm

Gross weight: 54KG

Net weight: 52KG

Option: I.V. Pole, Mattress



to order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663


2-Function up/down manual Stainless steel stretcher BD207-2

Features:

The emergency bed is used for carrying patients and the wounded people for the hospitals.

Material: high quality stainless steel and polishing .

With two function controlled by hand cranks to adjust the height and back section, which could make the operation more relaxing.



Specifications:

Standard Configuration	Four 5" silent castors 1 pair of stainless steel side rails 3 sets of cranks 1 piece of shelf 1 piece of Oxygen cart holder 2 Urine hooks. ABS foldable crank handles Switch of the side rails
Measurement	Size: 1950*700*450~680mm Gross weight: 54KG Net weight: 50KG Backrest tilting $\geq 70^\circ$ Hi-Low :450-680mm
Optional Accessories	I.V. Pole, Mattress

to order

All specifications, pictures are subject to change without earlier notification.
www.wanrooemed.com TEL: 0086-512-58810860 FAX: 086-512-58810663