



# Family Rehabilitation Nursing

## Rehabilitation Series

- Manual wheelchair/
- Electric wheelchair/
- Rollator/walker
- commode chair

See more happiness

FDA approved wheelchair with top quality  
Widely sold in USA, Europe, Japan and Korea

## Steel Wheelchair



SMW01 >>

Powder coated (Carbon) steel wheelchair, plaid nylon upholstery, fixed full armrest, 7" front rubber solid tire, 24" rear spoke wheel with pneumatic tire, fixed footrest with plastic Footplate.



SMW03 >>

Powder coated (Carbon) steel wheelchair, black nylon upholstery, Flip back and removable desk armrest, 8"\*1" front rubber solid tire, adjustable in 3 position , 24" rear mag wheel with PU tire, with height adjustable footrest and plastic footplate.

MODEL NO.				 (KG)	 DIM	 NW   GW	
SMW01	650	480	900	115	1020×260×920	17	20
SMW03	690	520	920	115	830×320×925	19	21



# Steel Wheelchair



**SMW04** >>

Powder coated (Carbon) steel wheelchair, black nylon upholstery(flame retardant), removable desk armrest, 8”\*1” front rubber solid tire, 24” rear quick release spoke wheel with pneumatic tire, E&J footrest with plastic footplate, folding back to reduce the box size.



**SMW20** >>

Chromed reclining wheelchair, raised headrest, black vinyl upholstery(flame retardant), Removable full armrest,200\*50mm front PU solid tyre, 24” rear spoke wheel with pneumatic tire and steel handrim , With elevating legrest , hard seat,commode and anti-tippers.It can










**SMW09** >>








Chromed Steel commode wheelchair, PU hard seat with commode, PU back, fixed armrest and footrest with heel loop, 24” rear spoke pneumatic wheel and 8” front caster ,optional Tray available.



**SMW30** >>

Chromed heavy duty steel wheelchair , double cross bar , black PU seat and back , flip up armrest and swing away footrest , 24” heavy duty rear mag wheel with plastic handrim, height adjustable swing away footrest with plastic footplate, 20”/22”/24” seat width optional.

MODEL NO.				 (KG)	 DIM	 NW	 GW
SMW04	690	520	920	115	850×320×930	19	21
SMW09	640	530	895	115	1020×320×920	19.2	23

MODEL NO.				 (KG)	 DIM	 NW	 GW
SMW20	670	510	1240	115	900×365×910	31.5	33.5
SMW30	770	505	910	130	830×460×930	25.8	29



# Aluminum Wheelchair








AMW01 >>

Aluminum wheelchair, anodized frame, double crossbrace, canvas upholstery, 8” front PU caster, 24”rear spoke wheels with pneumatic tyres,plastic handrim, fixed armrest and fixed footrest with plastic footplates and legband, with loop brake



AMW02 >>

Aluminum light wheelchair, flip up short desk armrest, swing away footrest, 24” rear spoke wheel with pneumatic tire, folding back, nylon seat and back, plastic footplate with heel loop.

MODEL NO.				 (KG)	 DIM	 NW   GW	
AMW01	645	490	925	115	1100×290×940	13.7	16.8
AMW02	660	500	940	115	830×330×940	16	19.2








AMW04 >>

Aluminum wheelchair, hand brake, fixed aremrest and footrest with heel loop, 6” front caster and 22” rear spoke wheel with pneumatic tire.



AMW05 >>

Aluminum light weight wheelchair, fixed armrest and footrest, 24” rear spoke wheel with pneumatic tire, plastic footplate with heel loop

MODEL NO.				 (KG)	 DIM	 NW   GW	
AMW04	630	470	870	115	780×290×910	11.2	14
AMW05	635	460	890	115	950×320×920	12.8	15.5





### AMW06 >>

Lightweight aluminum wheelchair, nodized frame, fixed armrest and height adjustable footrest, colorful plaid nylon upholstery and seat belt, 6" front casters and 12" rear spoke wheel with pneumatic tires or solid tire, with folding back option



### AMW09 >>

Completely anodized aluminum wheelchair, PVC upholstery, 8" front rubber tyres, 24" spoke wheels with pneumatic tyres, plastic driving wheel, fixed footrest with plastic footplates, with tray and commode.

MODEL NO.							
				(KG)	DIM	NW	GW
AMW06	590	450	900	115	680×290×710	10	12
AMW09	640	510	930	115	1040×280×960	17.8	20



### AMW10X >>

Deluxe alum wheelchair, flip back and removable full armrest, armrest height adjustable, back height adjustable, folding back option, 7"PU front caster, 12" or 24" quick release rear alum spoke rear wheel with PU tire, with loop brake for easy maneauvebility.



### AMW10XX>>

Deluxe alum wheelchair, flip back and removable full armrest, armrest height adjustable, back height adjustable, folding back option, 7"PU front caster, 12" or 24" quick release rear alum spoke rear wheel with PU tire, with loop brake for easy maneauvebility (expanding brake).

MODEL NO.							
				(KG)	DIM	NW	GW
AMW10X	675	510/560	945	130	850×360×850	16.5	19.6
AMW10XX	675	510	945	130	850×360×850	16.7	19.8



# Aluminum Wheelchair









**KS-A01** >>

Aluminum wheelchair, fixed long armrest and footrest, 6" front caster, 22" rear wheel



**KS-A02** >>

Aluminum functional wheelchair, elevating and flip up armrest, detachable footrest, adjustable folding back. 7" front casters, 24" rear wheels, seat width: 44cm, (6" front caster, 16" 、 20" rear wheel can be attached; 7" front caster, 20" 、 22" 、 24" rear wheel can be attached

MODEL NO.				 (KG)	 DIM	 NW    GW	
KS-A01	660	450	850	115	930×320×900	11.5	14
KS-A02	630	485	885	115	930×310×930	14	17.1





**KS-A03** >>

Deluxe aluminum alloy wheelchair, detachable armrest and footrest, adjustable orthopedics footrest and high-back (all the lying), with anti-tipper & commode. 8" front caster, 24" rear wheel (6" front caster, 16" rear wheel can be attached)



**KS-A04** >>

Lightweight aluminum wheelchair, total weight is only 11kgs, anodized frame, fixed armrest and footrest, colorful plaid nylon upholstery and seat belt, 6" front casters and 20" rear spoke wheel with solid tire,

MODEL NO.				 (KG)	 DIM	 NW    GW	
KS-A03	690	495	1215	115	910×340×945	21	23
KS-A04	585	450	845	100	930×300×865	11	13.4



# Aluminum Wheelchair



**KS-A05** >>

Aluminum frame, fixed armrest and footrest, folding back, aluminum front fork, With handbrake. Height adjustable footrest, Pneumatic tyre, 6” front wheel, 22” aluminum rear wheel,folding to save space, with padded seat and back.



**STC01** >> Steel Transport Wheelchair

Steel transport wheelchair, fixed full armrest, swing away footrest with plastic footplate, black nylon upholstery, all 8”\*1” solid tire, with folding back.

MODEL NO.				 (KG)	 DIM	 NW	 GW
KS-A05	675	440	905	100	940×355×710	14.3	17.1
STC01	570	490	970	100	625×260×810	12	14



**ATC01** >>

Aluminum lightweight transport chair, fixed full armrest, swing away footrest with plastic footplate, black nylon upholstery, all 8”\*1” rubber solid tire, with folding back.



**ATC02** >>

Aluminum lightweight transport chair, fixed full armrest, swing away footrest with plastic footplate, black nylon upholstery, 8”\*1” rubber solid tire, 12” rear PU tyre, with folding back option and loop brake.

MODEL NO.				 (KG)	 DIM	 NW	 GW
ATC01	570	480	1005	100	625×260×810	10.5	12
ATC02	570	480	1005	100	680×280×780	14	16



# Electric Wheelchair



KSE-01 >>

Half-reclining electric wheelchair, detachable seat, PU tyre, 33AH lead acid battery, adjustable width and height of armrest. Flip-up and folding footrest, magnetic brake. Convenient operation of electric control equipment. Buttoned seat belt. Shock absorption ways of Lever type.



KSE-02 >>

Alum. Frame with overall design of battery case, folding without removing battery, Flip-up armrest and swing away footrest, 16” rear PU tire and 200\*50 PU caster, magnetic brake, 20AH rubber battery, PG controller give you convenient and safety operation .

MODEL NO.				 (KG)	 DIM	 NW	 GW
KSE-01	620	520	1080	136	830×660×950	80	90
KSE-02	590	510	950	100	970×440×940	45.6	63



KSE-03 >>

High strength alum. folding frame, buckle type design of battery case. Easy to carry and store, 16” rear PU tire and 200\*50 front PU caster, magnetic brake, 20AH rubber battery, adjustable armrest height and width. PG controller give you convenient and safety operation .



KSE-04 >>

Alum. foldable frame, flip up armrest and swing away footrest, fantastic 16” rear wheel with build in driving motor , 8” front PU caster , magnetic brake, 20AH battery support 20miles driving , padded seat and safety belt, light weight and easy operation .

MODEL NO.				 (KG)	 DIM	 NW	 GW
KSE-03	680	480	905	100	970×440×940	43	48
KSE-04	690	500	935	100	725×430×945	43	49



# Electric Wheelchair



KSE-05 >>

Alum.foldable frame , flip up armrest and swing away footrest, special design longer backrest could lying down 180 degree, magnetic brake, 24” rear wheel with build-in driving motor , hand rim easy for manual using, 20AH battery support 20miles driving.type.Alum.foldable frame , flip up armrest and swing



KSE-07 >>

Light weight alum.frame foldable, 16” solid rear wheel and 7” PU caster, 20AH battery and magnetic brake easy operated controller indicating battery working station










KSE-06 >>








High strength alum.frame, light weight and foldable,16” solid rear wheel and 10” PU caster, magnetic brake, 20AH battery support 20miles driving, easy operated controller indicating battery working station



KSE-08 >>

Light weight alum. frame foldable, height adjustable and flip up armrest, swing away footrest , 16” pneumatic rear wheel and 7” PU caster, magnetic brake, 20 AH battery and 200W driving motor support manual / power model switching

MODEL NO.				 (KG)	 DIM	 NW	 GW
KSE-05	730	500	870	100	880×550×890	44.5	49
KSE-06	610	550	935	75	790×560×950	39	45

MODEL NO.				 (KG)	 DIM	 NW	 GW
KSE-07	590	495	950	100	790×560×930	37.5	42
KSE-08	610	540	925	75	790×430×795	38	42



# Rehabilitation series



Item No.:SR01  
Description: Wheels steel Rollator  
Features:  
Loop brake, PU seat, all 200\*50mm casters, handlebar height adjustable, with wire basket and tray. Weight Capacity: 250 lbs



Item No.:AR02  
Description: Four wheeled Aluminum Rollator  
Features:  
Loop brake, large padded PVC seat, all 6”\*1” casters, curved padded back, handlebar height adjustable, with tote bag or wire basket under seat. Weight capacity: 300 lbs



Item No.:AR04  
Description: Four wheeled Aluminum Rollator  
Features:  
Loop brake, large padded PVC seat, all 200\*40mm casters, curved padded back, handlebar height adjustable, with wire basket, Weight capacity: 300lbs



Item No.:AW01  
Description: Aluminum Walker  
Aluminum oxidized folding walker, may make “N” migration and stationary migration, double supporting tubes. Armrest width: 54cm, Total height: 82cm



Item No.:AW02  
Description: Aluminum Walker  
Aluminum oxidized folding walker, may take “N” migration ,double supporting tube, with two caster . Weight capacity:250lbs



Item No.:CM01  
Description: powder coated commode  
Powder coated frame, detachable, easy to carry. Weight Capacity: 250 lbs



Item No.:CM02  
Description: Commode Chair  
Folding commode chair, carbon steel frame with chrome finish, 6’ wheels can be attached, Weight Capacity: 375 lbs



Item No.:CM03  
Description: Commode Chair  
Detachable commode chair, carbon steel frame with chrome finish 6’ wheels, detach-able armrest and legrest, detachable soft seat. Weight Capacity: 250 lbs