



HBHKN-2000

INFANT RADIANT WARMER



Stock code [300314]

JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.



Brief Introduction

- Three control Modes: Pre-warm Mode, Manual Mode and Baby Mode control by micro-computer;
- Set Temperature and skin temperature can be displayed separately;
- Main body is made of stainless steel;
- Double sides phototherapy radiant box;
- Acrylic glass bassinet;
- APGAR timer;
- Multiple failure alarm indications;
- RS232 connector for data output;
- Disposable skin temperature sensor and SpO₂ are optional.

Standard Configuration

Warmer module, temperature controller, bassinet, phototherapy radiant box, and the stand.

Specification

Power requirement: AC220V-230V/50Hz or AC110-120V/50-60Hz, 850VA

Control Mode: Pre-warm Mode, Manual Mode and Baby Mode are controlled by micro-computer

Baby Mode temperature control range: 34.5-37.5°C

Skin temperature display range: 5°C-65°C

Accuracy of skin temperature sensor: ±0.2°C

Temperature uniformity of mattress: ≤2°C

Dominant wavelength of blue light: 420-470nm

Timing range of blue light: 0~9999hrs 59min

APGAR timer: audible and visible alarm when the device runs to 50" ~1' , 4' 50" ~5' , 9' 50" ~10'

Failure alarm: over temperature alarm, deviation alarm, sensor failure alarm, power failure alarm, setting alarm, checking alarm and so on.

Average value of general bilirubin radiation on the mattress effective range:

≥250uW/cm² (the upper phototherapy equipment)

≥600uW/cm² (the nether phototherapy equipment)

Bilirubin uniformity of radiating on the mattress effective range: >0.4

Environmental Requirement

Operating range: 18°C~30°C

Ambient air velocity rate: <0.3m/s

Product Packaging

Package: The whole unit is packed into one carton

Dimension: 940mm x 570mm x 1080mm

Gross Weight: 57kg



JIANGSU ROOE MEDICAL TECHNOLOGY CO.,LTD.

Add: Room1311,1313, Guotai New World Plaza, No19, East Renmin Road,Zhangjiagang, Jiangsu, China.

Tel:+86 512 58810663 Whatsapp/Mobile:+ 86-18962287398;15962399164

Email:info@wanrooemed.com Website: www.rooemed.com www.wanrooemed.com